# PHILIPS

## HF 305

# PHILIPS

**English****Page 4**

- Keep page 3 open when reading these operating instructions.

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
## Important

- Read these instructions for use carefully in order to obtain accurate measuring results.
- This blood pressure meter is only suitable for measuring the blood pressure of adults.
- This blood pressure meter measures blood pressure on the basis of the pulse rate. If you have a very weak or irregular pulse, it is possible that the blood pressure meter will not be able to determine the right blood pressure value.  
Diabetics and people who suffer from cardiac arrhythmia, vascular constriction or liver disorders must consult their doctor before they start measuring their own blood pressure, since their blood pressure values might deviate from the normal values. This also holds for pregnant women.
- If in doubt, consult your doctor.
- Only a doctor is capable of interpreting your blood pressure values; never change your medication without consulting a doctor.
- The cuff is suitable for upper arms with a circumference of 22-32 cm.


## List of functions

- A** Systolic pressure in mmHg (= mm mercury pressure)
- B** Diastolic pressure in mmHg (= mm mercury pressure)
- C** Pulse rate (the number of heartbeats per minute)
- D** Release button for the lid
- E** On/off button
- F** Memory button for previously measured value
- G** Protective lid with brief directions for use
- H** Hose connection
- I** Arm cuff
- J** Connector

## Placing the batteries


Place the batteries as shown in the figure on the left. Make sure that the + and - poles are pointing in the right direction. Use 4 alkaline, Philips Powerlife or rechargeable batteries, size AA (LR06).

Remove the protective foil from the display.


## Placing the cuff

Make sure that the blood circulation is not constricted by tight clothes, such as rolled-up shirt sleeves, as this could cause inaccurate readings.


Press the connector firmly into the hose connection of the appliance.


Put the cuff on a table as shown in the figure.

Put your left arm through the opening in the cuff. The rubber hose attached to the cuff should be on the inside of the arm, pointing in the direction of the hand.

Tighten the cuff around the upper arm in such a way that you can still slip two fingers between the cuff and the arm. Make sure that the bottom edge of the cuff is placed approx. 2-3 cm (1 inch) above the elbow


If the cuff has been fastened too tightly or not tightly enough, there is a possibility that the measurement will be inaccurate.

It is important that you sit down during the measurement and rest your forearm on a table. Keep the arm still during the measurement.

### Conditions for accurate blood pressure measurement

- Make sure you are comfortably seated and rest your forearm on a table in such a way that the upper arm is at about the same height as your heart.
- Sit down and relax for 10 minutes before you use the blood pressure meter for the first time. Excitement or depression as a result of emotional stress have a clear effect on blood pressure. Blood pressure will be higher (or lower) than normal and the pulse will usually also be faster.
- Make sure the measurement is not affected by tight clothes causing constricted blood vessels.
- Do not eat, drink or smoke before you measure your blood pressure. Alcoholic beverages in particular will influence the measurement.
- Carry out the measurement every day at the same time, since your blood pressure fluctuates in the course of the day.
- Do not attach too much importance to the results of a single measurement. Carrying out measurements regularly over a longer period of time and recording the results will give you a better insight into your blood pressure than a single measurement.
- If the blood pressure meter registers an abnormal situation, the error message 'Err' will appear on the display. For more details, see section 'Messages on the display'.


### Carrying out the measurement

Press the START button. All symbols will appear on the display for a second.

After the '0' has been displayed for a few seconds, the cuff will inflate automatically.


**Note:** The cuff should be fastened in such a way that the blue mark on the cuff is positioned over the main arteries in the upper arm.

When the correct blood pressure has been reached, pressurisation is stopped and the measurement starts (1). Remain seated in a relaxed way and do not move.

The pressure in the cuff will now be reduced. The 'measurement in progress' symbol ♥ will appear on the display (2). The measurement is completed as soon as the cuff is automatically depressurised.

Your systolic pressure, diastolic pressure and pulse will now appear on the display.

Press the START button again to switch the appliance off. If you forget to switch the appliance off, it will switch off automatically after one minute.

**Note:** wait approx. 5 minutes before you repeat a measurement.

## Memory


The values measured last are stored in the memory when the appliance is switched off.

The stored data appears on the display when the memory button is pressed. An 'M' is displayed in the bottom left-hand corner of the display.

The memory will be erased when you replace the batteries.

## Messages on the display

<b>Symbol</b>	<b>What to do</b>
	Measurement is in progress. Remain seated and do not move.
	The batteries are almost empty. Replace all batteries. Use alkaline, Philips Powerlife or rechargeable batteries, size AA (LR06).
	This message is displayed if you have moved during the measurement. Measure your blood pressure again and do not move during the measurement.
	The pressure in the cuff is not high enough. Check if the hose has been properly connected to the monitor and if the cuff has been placed around the arm correctly. Check whether there is any air escaping from the cuff.

## General information about blood pressure

Blood pressure is measured in mmHg (mm mercury pressure). Blood pressure consists of two values. Systolic pressure (SYS) is created by the contractions of the heart muscle which pumps the blood into the arteries. Diastolic pressure (DIA) is created when the heart relaxes and fills with blood again. The World Health Organisation (WHO) has laid down the following standard values with regard to the blood pressure of adults:

### High blood pressure = risk factor

Systolic pressure (SYS) is too high when it exceeds 160 mmHg and diastolic pressure (DIA) is too high when it exceeds 95 mmHg. We speak of 'high blood pressure' when one of the two pressures exceeds the value mentioned above.

### Borderline values

Systolic pressure: 140-160 mmHg

Diastolic pressure: 90-95 mmHg

### Standard values

Systolic pressure: max. 140 mmHg

Diastolic pressure: max. 90 mmHg

## Causes and consequences of high blood pressure


Causes of high blood pressure are, among other things, obesity, a high cholesterol level, smoking, excessive alcohol consumption, emotions, excessive salt consumption, a lack of exercise, hereditary inclination towards high blood pressure and certain diseases such as kidney and metabolic disorders.

In the long term, high blood pressure can lead to reduced elasticity of the walls of arteries (arteriosclerosis), which is one of the major causes of cardiac infarcts, cerebral haemorrhages and cerebral infarcts as well as kidney disorders.

## Blood pressure fluctuation

The blood pressure of healthy people, but especially that of people with high blood pressure, fluctuates greatly. Blood pressure fluctuates in the course of the day and in the course of the month. The season and the temperature also have an impact on blood pressure. Generally, a person's blood pressure is lower when the weather is hot and higher when it is cold.

To give an example of the fluctuations in blood pressure that occur during the day, a graph showing the blood pressure of a 35-year-old man has been included. The blood pressure fluctuations may amount to 30 to 50 mmHg as a result of changing circumstances during the day. Normally, blood pressure will go up when a person is working or exercising and will drop to its lowest level during the hours of sleep. So there is no need to be alarmed if one measurement shows your blood pressure to be higher than usual.

## Why measure blood pressure at home?

An incidental blood pressure measurement carried out by a doctor is not a true representation of your actual blood pressure. It is possible that a person's blood pressure is higher when it is measured by a doctor than it is when it is measured at home, because people tend to be somewhat tense during a visit to the doctor. This is why the pressure measured by doctors is often called 'white coat' blood pressure. By measuring your blood pressure every morning before breakfast you can determine your basic blood pressure. By always carrying out the measurement at the same time of day, external influences on blood pressure can be excluded as much as possible.

If you want to gain insight into the daily fluctuations of your blood pressure, measure your blood pressure at several set points during the day. The risk of illness due to high blood pressure depends on all pressures that occur during the day, not only on the blood pressure you happen to have in the morning.

Always write down the date and the time along with the blood pressure values measured. Keep in mind that only a doctor is capable of analysing your blood pressure. On no account make a diagnosis yourself and NEVER change the dose of medication prescribed by a doctor.

## Maintenance

- This digital blood pressure meter consists of very fragile precision parts. Be careful to avoid extreme temperatures, humidity, sunlight, shocks and dust.
- Only clean the appliance with a soft, moist cloth and some washing-up liquid. Carefully remove stains on the cuff with a moist cloth and some soap.
- Do not bend the cuff too far backwards and make sure that the cuff has not been rolled up too tightly when you put it away.
- Remove the batteries from the appliance if the blood pressure meter will not be used for quite some time.
- Do not throw the batteries away but hand them in at an officially assigned collection point.

## Technical data Philips HF 305

<b>Measuring method</b>	<b>Oscillometric</b>
Measuring range	Pressure: 20 - 280 mmHg Pulse: 40 - 200 beats/minute
Measuring accuracy	Pressure: $\pm 3$ mmHg Pulse: $\pm 5\%$
Batteries	4 x 1.5V alkaline or Philips Powerlife batteries (AA/LR06)
Battery life	> 300 measurements
Operating temperature	+10°C to +40°C.
Storage temperature	-10°C to +60°C.
Air humidity	85% or less
Weight	approx. 320 g without batteries

*This appliance meets the requirements of MDD 93/42/EEC for medical appliances*

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