PHILIPS

HF 320

PHILIPS

English**Page 4**

- Keep page 3 and open when reading these operating instructions.

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

Important

- **Read these instructions for use carefully in order to obtain accurate measuring results.**
- **This blood pressure meter is only suitable for measuring the blood pressure of adults.**
- **This blood pressure meter measures blood pressure on the basis of the pulse rate. If you have a very weak or irregular pulse, it is possible that the blood pressure meter will not be able to determine the right blood pressure value.**

Diabetics and people who suffer from cardiac arrhythmia, vascular constriction or liver disorders must consult their doctor before they start measuring their own blood pressure, since their blood pressure values might deviate from the normal values. This also applies to pregnant women.

- **If in doubt, consult your doctor.**
- **Only a doctor is capable of interpreting your blood pressure values; never change your medication without consulting a doctor.**

List of functions



(See picture on the inside of the sleeve)

- A** Systolic pressure in mmHg (= mm mercury pressure)
- B** Diastolic pressure in mmHg (= mm mercury pressure)
- C** Pulse rate (the number of heartbeats per minute)
- D** Release button for the lid
- E** On/off button
- F** Memory button for seven previously measured values
- G** Protective lid with brief directions for use
- H** Wristband

Before the first measurement

- Remove the protective foil from the display.
- Place the two AAA (LR03) batteries in the battery compartment. Make sure that the + and - poles are pointing in the right direction.

Please note: Insert the lower battery + pole first.



Placing the wristband

- Put the wristband on the left wrist with the control panel on the inside of the wrist. Make sure that there is no piece of clothing between the wrist and the wristband.

The wristband must be placed approx. 1 cm (0.5 inch) away from the palm of the hand

- Tighten the wristband firmly around the wrist; make sure that there is no room left between the wrist and the wristband.

Note: the coloured marks on the wristband indicate a wrist circumference of 13.5-19.5 cm (5.25-7.75 inches). If the marks do not overlap, the circumference of your wrist is too big to allow accurate functioning of this blood pressure meter.

Conditions for accurate blood pressure measurement

- Rest your elbow on a table and hold up your forearm in such a way that the wristband is at about the same height as your heart. Support your forearm.
- Sit down and relax for 10 minutes before you use the blood pressure meter for the first time. Breathe deeply in and out five times before you start measuring.

- You have to wait at least 3 minutes between 2 measurements to avoid incorrect values caused by temporary congestion of the arteries.
- Carry out the measurement every day at the same time, since your blood pressure fluctuates in the course of the day.
- Do not eat, smoke or drink before you measure your blood pressure. Alcoholic beverages in particular will influence the measurement.
- Do not attach too much importance to the results of a single measurement. Carrying out measurements regularly over a longer period of time and recording the results will give you a better insight into your blood pressure than a single measurement.



Carrying out the measurement

- Sit down comfortably, put the wristband around your wrist and hold it at the same height as your heart.

Note: if you expect your systolic pressure to exceed 150, read the instructions in chapter "Systolic pressure higher than 150 mmHg".

- Press the START button. All symbols will appear on the display for a moment.
- The wristband will inflate automatically and the display will show the increasing pressure.

When the correct pressure has been

reached in the wristband, pressurisation is stopped and the measurement starts.

Note: *the blood pressure meter will automatically increase the pressure in the wristband if the pressure in the wristband turns out to be too low during the measurement. Remain seated in a relaxed way.*

- During the measurement the pressure in the wristband will be slowly reduced and the ♥ symbol will appear to indicate that the measurement is in progress.
- When the measurement has been completed, the cuff is automatically depressurised. Your systolic pressure (SYS), diastolic pressure (DIA) and pulse (PUL) will now appear on the display.

You can switch the appliance off by pressing the START button again. If you forget to switch the appliance off, it will switch off automatically after three minutes.

Memory

The appliance automatically stores previously measured values in the memory. Up to seven measurements can be stored. You can retrieve this data in the following way:

- Press the MEMORY button; the data of the most recently stored measurement is displayed.

- Every time you press the MEMORY button, the data of the measurement previous to the one currently shown will appear on the display.
- The memory will be erased when you replace the batteries.

Systolic pressure higher than 150 mmHg

If you switch the appliance on in the normal way, the wristband will be inflated to 180 mmHg. If your systolic pressure is higher than 150 mmHg, the appliance will have to interrupt the measurement in progress in order to inflate the wristband to the correct pressure. This will cause your wrist to be squeezed tight for quite some time. You can avoid this inconvenience by following the instructions below, which will make the wristband inflate to the correct value straight away.

- Put the wristband around your wrist.
- Sit down comfortably and hold the wristband at the same height as your heart.
- Press the start button and keep it pressed in until the value shown on the display is 40 mmHg higher than your expected systolic pressure.

Example: if you expect your systolic pressure to be 160 mmHg, keep the START button pressed in until the value shown on the display is 200.

Release the START button.

- The measurement will now take place fully automatically as described in chapter “Carrying out the measurement”.

Messages on the display

<i>Symbol</i>	<i>What to do</i>
	Measurement is in progress. Remain seated and do not move.
<i>Err</i>	This message is displayed if you have moved during the measurement. Measure your blood pressure again and do not move during the measurement.
<i>Err</i> <i>CUF</i>	The wristband has not been put around the wrist correctly. Take it off and put it back on again and carry out the measurement once more.
	The batteries are almost empty. Replace both batteries. Use alkaline, Philips Powerlife or rechargeable batteries, size AAA (LR03).

General information about blood pressure

Blood pressure is measured in mmHg (mm mercury pressure). Blood pressure consists of two values. Systolic pressure (SYS) is created by the contractions of the heart muscle which pumps the blood into the arteries. Diastolic pressure (DIA) is created when the heart relaxes and fills with blood again. The World Health Organisation (WHO) has laid down the following standard values with regard to the blood pressure of adults:

High blood pressure = risk factor

Systolic pressure (SYS) is too high when it exceeds 160 mmHg and diastolic pressure (DIA) is too high when it exceeds 95 mmHg.

We speak of 'high blood pressure' when one of the two pressures exceeds the value mentioned above.

Borderline values

Systolic pressure: 140-160 mmHg

Diastolic pressure: 90-95 mmHg

Standard values

Systolic pressure: max. 140 mmHg

Diastolic pressure: max. 90 mmHg

Causes and consequences of high blood pressure

Causes of high blood pressure are, among other things, obesity, a high cholesterol level, smoking, excessive alcohol consumption, emotions, excessive salt consumption, a lack of exercise, hereditary inclination towards high blood pressure and certain diseases such as kidney and metabolic disorders.


In the long term, high blood pressure can lead to reduced elasticity of the walls of arteries (arteriosclerosis), which is one of the major causes of cardiac infarcts, cerebral haemorrhages and cerebral infarcts as well as kidney disorders.

Blood pressure fluctuation

The blood pressure of healthy people, but especially that of people with high blood pressure, fluctuates greatly. Blood pressure fluctuates in the course of the day and in the course of the month. The season and the temperature also have an impact on blood pressure. Generally, a person's blood pressure is lower when the weather is hot and higher when it is cold.

To give an example of the fluctuations in blood pressure that occur during the day, a graph showing the blood pressure of a 35-year-old man has been included.

The blood pressure fluctuations may amount to 30 to 50 mmHg as a result of changing circumstances during the day. Normally, blood pressure will go up when a person is working or exercising and will drop to its lowest level during the hours of sleep. So there is no need to be alarmed if one

measurement shows your blood pressure to be higher than usual.

Why measure blood pressure at home?

An incidental blood pressure measurement carried out by a doctor is not a true representation of your actual blood pressure. It is possible that a person's blood pressure is higher when it is measured by a doctor than it is when it is measured at home, because people tend to be somewhat tense during a visit to the doctor. This is why the pressure measured by doctors is often called 'white coat' blood pressure.

By measuring your blood pressure every morning before breakfast you can determine your basic blood pressure. By always carrying out the measurement at the same time of day, external influences on blood pressure can be excluded as much as possible.

If you want to gain insight into the daily fluctuations of your blood pressure, measure your blood pressure at several set points during the day. The risk of illness due to high blood pressure depends on all pressures that occur during the day, not only on the blood pressure you happen to have in the morning.

Always write down the date and the time along with the blood pressure values measured. Keep in mind that only a doctor is capable of analysing your blood pressure values. On no account make a diagnosis yourself and NEVER change the dose of medication prescribed by a doctor.

Difference between measurement at upper arm and wrist

- The blood pressure measured with a watch-type blood pressure monitor may differ from the values measured on the upper arm. Measurement with a wrist blood pressure monitor may show a difference of up to 10 mmHg for a healthy adult. The blood vessels in the upper arm are connected to those in the wrist and, as a result, blood pressure fluctuations also occur in the wrist. Therefore, measuring blood pressure at the wrist is a good way of monitoring your blood pressure.

In those who suffer from hypertension, arteriosclerosis, diabetes, kidney failure and other disorders of the circulation system, the difference may be greater than 10 mmHg. In that case measurement on the upper arm is to be preferred.

In case of doubt, consult your doctor about the best method in your case or make a comparison between a measurement on the upper arm and one on the wrist.

Maintenance

- This digital blood pressure meter consists of very fragile precision parts. Be careful to avoid extreme temperatures, humidity, sunlight, shocks and dust.
- Only clean the appliance with a soft, moist cloth and some washing-up liquid.
Carefully remove stains on the wristband with a moist cloth and some soap.
- Do not bend the wristband too far backwards and make

sure that the wristband has not been rolled up too tightly when you put it away.

- Remove the batteries from the appliance if the blood pressure meter will not be used for quite some time.
- Do not throw the batteries away but hand them in at an officially assigned collection point.

Technical data Philips HF 320

Measuring method	Oscillometric
Measuring range	Cuff pressure: 20 - 280 mmHg Pulse: 40 - 200 beats/minute
Measuring accuracy	Cuff pressure: ± 3 mmHg Pulse: $\pm 5\%$
Memory	7 measurements
Batteries	2 x 1.5V alkaline or Philips Powerlife batteries (AAA/LR03)
Wristband	Wrist circumference of 13.5-19.5 cm. (5.25 - 7.75 inches)
Operating temperature	+10°C to +40°C.
Storage temperature	-20°C to +70°C.
Air humidity	85% or less
Weight	approx. 140 g without batteries

This appliance meets the requirements of MDD 93/42/EEC for medical appliances.

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