PHILIPS

blood pressure monitor
HF321, 319

ENGLISH 4

NORSK 17

SVENSKA 30

SUOMI 43


DANSK 55

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Important


- ▶ Read these directions for use carefully in order to obtain accurate measuring results.
- ▶ This blood pressure monitor is only suitable for measuring the blood pressure of adults.
- ▶ This blood pressure monitor measures blood pressure on the basis of the pulse rate. If you have a very weak or irregular pulse, the blood pressure monitor may not be able to determine the right blood pressure value.
- ▶ Diabetics and people who suffer from cardiac arrhythmia, vascular constriction or liver disorders must consult their doctor before they start measuring their own blood pressure, since their blood pressure values might deviate from normal values. Pregnant women should also consult their doctor first.
- ▶ If in doubt, consult a doctor.
- ▶ Only a doctor is capable of interpreting your blood pressure values. Never change your medication without consulting a doctor.
- ▶ Your blood pressure may fluctuate during the day and it may even vary within minutes. This means that successive measurements may show different values.
- ▶ If the blood pressure monitor has been stored at a very low temperature (close to freezing point), put it in a warm place for at least one hour before using it.

Before the first measurement

1 Remove the protective foil from the display.

2 Open the battery compartment and place 2 AAA (LR03) batteries in the compartment.

3 Make sure the + and - poles are pointing in the right direction.

Push the + pole in first.


Measuring the blood pressure

Conditions for accurate blood pressure measurement

- ▶ Sit down and relax for 10 minutes before you use the blood pressure monitor. Breathe deeply in and out five times before you start measuring.
- ▶ Wait at least 3 minutes between 2 measurements to avoid incorrect readings due to the temporary congestion of the arteries.
- ▶ Perform the measurement every day at the same time, since your blood pressure fluctuates in the course of the day.
- ▶ Do not eat, smoke or drink before you measure your blood pressure. Alcoholic beverages in particular will influence the measurement.
- ▶ Do not attach too much importance to the results of a single measurement. Performing

measurements regularly over a longer period of time and recording the results will give you a better insight into your blood pressure than a single measurement.


Performing the measurement


- 1 Put the wristband on the left wrist with the control panel on the inside of the wrist. Make sure there is no piece of clothing between the wrist and the wristband.

The wristband must be placed approx. 1 cm (0.5 inch) away from the palm of the hand.

- 1 Tighten the wristband firmly around the wrist. Make sure there is no room left between the wrist and the wristband.

- The coloured marks on the wristband indicate a wrist circumference of 13.5 - 22 cm (5.3 - 8.3 inches). If the marks do not overlap, the circumference of your wrist is too large or too small to allow this blood pressure monitor to function accurately.

- 2 Sit down comfortably and hold your wrist at the same height as your heart.

If you expect your systolic pressure to exceed 150 mmHg, read the instructions under the heading 'Systolic pressure higher than 150 mmHg'.

3 Press the START button.

All symbols appear on the display for a moment. The wristband inflates automatically and the display shows that the pressure is increasing. When the correct pressure in the wristband has been reached, inflation of the wristband stops and the measurement starts.

The blood pressure monitor will automatically increase the pressure in the wristband if the pressure in the wristband turns out to be too low during the measurement. Remain seated in a relaxed way.

► The pressure in the wristband will slowly decrease during the measurement.

This is shown by the counting down of the pressure on the display.

Note: do not talk, cough, laugh or move your arm or hand muscles during the measurement.

► When the measurement has been completed, the cuff is automatically deflated. Your systolic pressure (SYS), diastolic pressure (DIA) and pulse (PUL) are now shown on the display.

Systolic pressure higher than 150 mmHg

If you switch the appliance on in the normal way, the wristband will be inflated to 180 mmHg.

If your systolic pressure is higher than 150 mmHg, the appliance will have to interrupt the measurement in progress in order to inflate the wristband to the correct pressure. This will cause your wrist to be squeezed tight for quite some time. You can avoid this inconvenience by inflating the wristband to the correct value straight away.

1 To inflate the wristband to the correct value: keep the **START** button pressed in until the value shown on the display is 40 mmHg higher than your expected systolic pressure.

► *Example: if you expect your systolic pressure to be 160 mmHg, keep the Start button pressed in until the value shown on the display is 200.*

Memory function

The blood pressure monitor automatically stores previously measured values in its memory.

Up to 60 measurements can be stored.

1 Retrieve data from the memory by pressing the **MEMORY** button.

The data of the last measurement are displayed.

2 Press the **MEMORY** button again and the data of the measurement previous to the one currently shown will appear on the display.

no 138
60 80
□ ▾ 77

no 123
59 76
□ ▾ 65

If you replace the batteries while the appliance is switched on, all data stored in the memory will be erased.

Cleaning and maintenance

- ▶ This digital blood pressure monitor consists of very fragile precision parts. Be careful to avoid extreme temperatures, humidity, sunlight, shocks and dust.
- ▶ Only clean the appliance with a soft, moist cloth and some washing-up liquid.
- ▶ Carefully remove stains on the wristband with a moist cloth or some soap.
- ▶ Do not bend the wristband too far backwards and do not roll the wristband up too tightly when you put it away.
- ▶ Remove the batteries from the appliance if the blood pressure monitor is not going to be used for quite some time.

Environment

- 1 Do not throw the batteries away with the normal household waste. Hand them in at an officially assigned collection point.

Technical data Philips HF321, HF319

Measuring method: Oscillometric

Measuring range: Pressure: 20 - 280 mmHg

Pulse: 40 - 180 beats/minute

Measuring accuracy: Pressure; +/- 3 mmHg

Pulse: +/- 5%

Memory: 60 measurements

Batteries: 2 x 1.5 alkaline or Philips Powerlife batteries (AAA/LR03)

Wristband: Wrist circumference of 13.5 - 22 cm (5.3 - 8.3 inches)

Operating temperature: +10° to +40°

Storage temperature: -5° to +60°

Air humidity: 95% or less

Weight: approx. 130 grams without batteries

Automatic shut-off: after 30 seconds.

- ▶ *This appliance meets the requirements of MDD 93/42/EEC for medical appliances.*

If you need information or if you have a problem, please visit the Philips Web site at www.philips.com or contact the Philips Customer Care Centre in your country (you will find its phone number in the worldwide guarantee leaflet). If there is no Customer Care Centre in your country, turn to your local Philips dealer or

contact the Service Department of Philips Domestic Appliances and Personal Care BV.

General information about blood pressure

Blood pressure is measured in mmHg (mm mercury pressure). Blood pressure consists of two values. Systolic pressure (SYS) is created by the contractions of the heart muscle, which pumps the blood into the arteries. Diastolic pressure (DIA) is created when the heart relaxes and fills with blood again. The World Health Organisation (WHO) has laid down the following standard values with regard to the blood pressure of adults

High blood pressure is a risk factor

Systolic pressure (SYS) is too high when it exceeds 160 mmHg and diastolic pressure (DIA) is too high when it exceeds 95 mmHg.

We speak of 'high blood pressure' when one of the two pressures exceeds the value mentioned above.

Borderline values

Systolic pressure: 140-160 mmHg

Diastolic pressure: 90-95 mmHg

Standard values

Systolic pressure: max. 140 mmHg

Diastolic pressure: max. 90 mmHg

Causes and consequences of high blood pressure


Causes of high blood pressure are, among other things, obesity, a high cholesterol level, smoking, excessive alcohol consumption, emotions, excessive salt consumption, a lack of exercise, hereditary inclination towards high blood pressure and certain diseases such as kidney and metabolic disorders.

In the long term, high blood pressure can lead to reduced elasticity of the walls of arteries (arteriosclerosis), which is one of the major causes of cardiac infarcts, cerebral haemorrhages and cerebral infarcts as well as kidney disorders.

Blood pressure fluctuation

The blood pressure of healthy people, but especially that of people with high blood pressure, fluctuates greatly. Blood pressure fluctuates in the course of the day and in the course of the month. The season and the temperature also have an impact on blood pressure. Generally, a person's blood pressure is lower when the weather is hot and higher when it is cold.

To give an example of the fluctuations in blood pressure that occur during the day, a graph showing the blood pressure of a 35-year-old man has been included.

The blood pressure fluctuations may amount to 30 to 50 mmHg as a result of changing circumstances during the day. Normally, blood pressure will go up when a person is working or exercising and will drop to its lowest level during the hours of sleep. So there is no need to be alarmed if one measurement shows your blood pressure to be higher than usual. Fluctuations of 30 mmHg may even occur within minutes.

Why measure blood pressure at home?

An incidental blood pressure measurement carried out by a doctor is not a true representation of your actual blood pressure. It is possible that a person's blood pressure is higher when it is measured by a doctor than it is when it is measured at home, because people tend to be somewhat tense during a visit to the doctor. This is why the pressure measured by doctors is often called 'white coat' blood pressure.

By measuring your blood pressure every morning before breakfast you can determine your basic blood pressure. By always carrying out the measurement at the same time of day, external influences on blood pressure can be excluded as much as possible.

If you want to gain insight into the daily fluctuations of your blood pressure, measure your blood pressure at several set points during the day. The risk of illness due to high blood pressure depends on all pressures that occur during the day, not only on the blood pressure you happen to have in the morning.

Always write down the date and the time along with the blood pressure values measured. Keep in mind that only a doctor is capable of analysing your blood pressure values. On no account make a diagnosis yourself and NEVER change the dose of medication prescribed by a doctor.

Difference between upper-arm and wrist measurement

The blood pressure measured with a watch-type blood pressure monitor may differ from the values measured on the upper arm. Measurement with a wrist blood pressure monitor may show a difference of up to 10 mmHg in a healthy adult.

In those who suffer from hypertension, arteriosclerosis, diabetes, kidney failure and other disorders of the circulation system, the difference may exceed 10 mmHg. In that case measurement on the upper arm is to be preferred.

In case of doubt, consult your doctor about the best method in your case or make a comparison between a measurement on the upper arm and one on the wrist.

Troubleshooting

Symbol	Reason	Action
E1	This message is displayed if you have moved during the measurement.	Measure your blood pressure again and do not move during the measurement. If this error message continues to appear, take the appliance back to your dealer for repair.
E1	E1 will also appear if the pressure is higher than 300 mmHg, in that case it will be followed by immediate auto-deflation.	Measure your blood pressure again and do not move during the measurement. If this error message continues to appear, take the appliance back to your dealer for repair.
E2	This message will appear if one or more components malfunction.	Take the appliance back to your dealer for repair.
EE	This message will appear when the pressure drops to a level lower than 40 mmHg (due to weak signals).	Measure your blood pressure again and do not move during the measurement.
Err Cuf	The wristband has not been put around the wrist correctly.	Put the cuff around the wrist correctly and try again.
Err Cuf	This message may also appear in case of a malfunctioning cuff or pump.	If this error message continues to appear, take the appliance back to your dealer for repair.
	The batteries are almost empty.	Replace both batteries. Use alkaline, Philips Powerlife or rechargeable batteries, size AAA (LR03).

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Personal Care B.V.