**Cellesse**  
*Cellulite massage system*

HP 5230

**PHILIPS**

# Cellesse

<i>English</i>	<i>Page 3</i>
<i>Français</i>	<i>Page 17</i>
<i>Deutsch</i>	<i>Seite 31</i>
<i>Italiano</i>	<i>Pagina 45</i>
<i>Nederlands</i>	<i>Pagina 59</i>

## *English*

<i>Welcome to Cellesse</i>	<i>Page 4</i>
<i>Cellulite - a feminine problem</i>	<i>Page 5</i>
<i>Cellulite - can treatment help?</i>	<i>Page 7</i>
<i>The Alternative Treatment Options</i>	
<i>Cellesse - your personal Cellulite Massage System</i>	<i>Page 8</i>
<i>Cellesse in Action</i>	<i>Page 8</i>
<i>Clinical Tests show results</i>	<i>Page 10</i>
<i>Read carefully before using Cellesse</i>	<i>Page 11</i>
<i>Instructions for use</i>	<i>Page 13</i>

# Welcome to Cellesse

*... and your new campaign against Cellulite*

*Cellesse is a special massage system designed to help you reduce or prevent the unsightly appearance of Cellulite.*

*Cellesse is the result of years of research in the Philips Personal Care Institute. We have investigated the causes of Cellulite, tested different types of treatment available and developed a solution which will help you get the results you want.*

*In the following pages we look, in some detail, at the causes of Cellulite and how Cellesse can help you with this particularly feminine problem.*

*Note: carefully read all pages of this booklet, especially the chapter "Read carefully before using Cellesse".*

# *Cellulite- a feminine problem*

*So-called Cellulite is the increased visibility of the lumpy fat deposits through the surface of the skin. A common description is 'orange peel skin' because of its unique pitted, bumpy appearance.*

*Cellulite is a cosmetic beauty problem which affects the appearance of certain parts of your body, notably the hips, thighs and buttocks. In more extreme cases, Cellulite can be uncomfortable, even painful when touched or squeezed.*

*Cellulite should not be confused with "cellulitis" which is a medical term linked to bacterial infection of the skin and associated inflammation of the tissue.*

*There are two main influences to the formation of Cellulite:*

## *1) Natural mattress structure in the female skin*

*Firstly, Cellulite is due to the specific, however natural mattress structure of the female skin. It is a fact of nature that women have, unlike men, vertical fat deposits in their skin (men's are more horizontal).*

*This mattress structure can be made especially visible when the skin is squeezed.*

*All women have the possibility of developing Cellulite because of this structure of the female fat deposits.*


### *Try the Squeeze Test*

*This is easy to identify. Just gently squeeze an area of your upper thigh with both hands and you will most probably see the mattress shaped fat layer in your skin.*


*Especially when the fat deposits are full due to weight gain or when the skin is naturally thin (or becomes thin through the normal ageing process) the lumpy appearance of the skin can become more visible.*

### *2) Fluid retention*

*This natural skin condition can however progress to a stage in which Cellulite becomes visible on the legs without squeezing.*

*Vertical fat deposits in the female skin*

*Horizontal fat deposit in the male skin*

*In most of these cases fluid retention is affecting the fatlayer. Women tend to retain more fluid than men, especially in their legs. This water retention also affects the fat tissue in the skin resulting in excess fluid around the fat cells. This adds to the puffiness of the mattress and make the Cellulite layer swell.*

*If this fluid retention is allowed to continue for months and years, the fat deposits thicken and harden, and eventually hard lumps are being formed. These lumps can then be seen through the skin without squeezing.*

*This chronic situation can be painful when squeezed or pinched, and is, unlike the light and medium forms of Cellulite, almost irreversible.*

#### *Try the Pinch Test*

*To check whether you have problems with fluidretention in your legs you can gently pinch the skin of your upper thigh and try to lift it away from the underlying layers.*

*If this hurts, especially when you then press it, it could be a first indication of water retention and additional swelling of the Cellulite layer.*

*Two simple tests to determine what stage of Cellulite you have.*

*Top left: The Pinch Test*

*Above: The Squeeze Test*

*Left: Cellulite*

# *Cellulite can treatment help?*

*Once established, Cellulite will not disappear or reduce of its own accord. On the contrary, it is likely to increase. The only solution is to make the fight against Cellulite part of your regular personal care routine - and begin treatment before the problem develops any further than it has today.*

*The good news is that Cellulite can be reduced, but it can never be completely eliminated and will return when neglected.*

*Your fight against Cellulite may have already taken many forms: diet, exercise, the use of special creams and gels, massage and even surgery. Now Cellese adds an effective, purpose-designed approach which gets right down to the root causes.*

## *The alternative treatment options*

*Diet and Cellulite. As described before Cellulite is not necessarily associated with being overweight only. Dieting tends to reduce fat which will help to reduce the swelling of the fat layer. However "feminine curves" like buttocks and thighs are especially resistant to losing fat.*

*Next to dieting, there is evidence that certain foods are linked to the formation of Cellulite, especially to the water retention: tea,*

*coffee, alcohol and refined foods such as white flour and sugar should be avoided as much as possible.*

*Exercise can also help in fighting Cellulite. Low intensity aerobic exercise such as walking or swimming will help boost the blood and lymphatic circulation. However, to be effective, exercise needs to be localised, targeting areas of the female body where Cellulite is concentrated and this is not easy to achieve.*

*Massage. Next to diet and exercise, removing the fluid from the Cellulite layer through massage can also help reduce the swelling and the bumpiness.*

*Intensive localised massage can mobilise the fluid in the fat layer and in this way help the body to remove it through the normal circulation and lymphatic drainage systems. In this way the swelling of the Cellulite layer can be reduced.*

*When this is done regularly it can also prevent further swelling and hardening of the fat layer. So removing excess fluid now from the Cellulite layer should prevent future swelling and hardening of this layer.*


# *Cellesse*

## *Your personal Cellulite Massage System*

*For this purpose Philips has developed Cellesse, a compact hand-held appliance designed specifically to combat and reduce Cellulite. Cellesse directly addresses the Cellulite layer and can be used in the privacy of your own home.*

*Cellesse gives intensive, deep massage of the Cellulite layer in the skin to mobilise the excess fluid which is then removed naturally by the body.*

*The result is a reduction in the swelling and the bumpiness of the Cellulite layer.*

# *Cellesse in Action*


*Cellesse has been developed especially for those parts of the body prone to Cellulite and for this has been equipped with an effective massage system: Vacumotion.*

*Vacumotion is a continuous lifting and rolling movement which gives a deep and intensive massage of the Cellulite layer, with the least of friction. It consists of two rollers and a middle vacuum section which are electronically connected to result in one fluid massage action.*


*The process resembles a wave-like movement of the skin as the Cellesse unit is moved over the affected area. The direction for use is always towards the heart, i.e. upwards along the leg and affected areas are covered in a series of adjacent, parallel movements.*

### *Sensor control*


*A special sensor detects proper contact with the skin, only when this proper contact is made does the Vacumotion begin. This ensures effective treatment at all times.*

Celesse starts automatically when placed on the skin. The front roller gently rolls the skin into the vacuum area.

The skin is then lifted even further by means of vacuum.

The rear roller massages the lower layers of the elevated skin and then gently rolls it off.

### Five Intensity settings

Five degrees of intensity are available. As you get to know Celesse better, you will be able to find the setting that suits you best. We recommend to start in the beginning with the lowest setting and gradually increase depending on your specific skin condition and preference. As the skin gets more used to the treatment it will accept a higher setting.

### Vacu-Stop

Vacu-Stop is a simple and safe way of stopping the treatment at any time.

The stop switch is easily accessible from either side of the appliance.

### Designing an Effective Treatment Programme

The recommended frequency of use for Celesse is 10-15 minutes on each leg, three times a week. In consumer tests we have found that through this type of regular use in already one month a significant smoothening of the Cellulite layer inside the skin can be achieved. Of course results will differ from woman to woman, depending on the extent and severity of the Cellulite and your lifestyle.


We found that in most cases, results can be achieved, however regular use of Celesse is required as fluid retention will always return and put the fatlayer under pressure again.

# Clinical Tests show results

Clinical tests of the Cellesse system have been conducted under controlled conditions to assess the effectiveness of the treatment. Ultrasound scans of the inner skin were used to measure the effects (see pictures).

The results have confirmed that Cellesse can achieve a smoothening of the Cellulite layer inside the skin: on average 33% smoothening of the Cellulite layer after 1 month of use. Up to almost 50% may be achieved after two months based on the recommended frequency of use (3x per week; 15 minutes per leg).

Ultrasound scans: before and after Cellesse


Results of Cellesse using ultrasound scanning  
Average smoothening of the Cellulite layer inside the skin.

After 1 month: 33%

After 2 months: 49%

After 3 months: 55%

Graph showing the stages of reduction

# *Read carefully before using Cellesse*

*If treatment is discontinued, the swelling of the Cellulite layer will gradually return as fluidretention will continue in the female body. As an indication for this we found in our study that one month after treatment the Cellulite layer became less smooth again. However still 39% of smoothening was retained and two months after treatment 25% smoothening was still visible on the ultrasound scans.*

*This implies that regular usage of Cellesse will have the largest impact especially early in the treatment cycle.*

*Therefore, we recommend intensive usage for the first 2-3 months and, after that, less intensive usage to keep fluid retention under control: for example, 15 minutes every week or 5 minutes 3 times per week, depending on your preference and skin condition.*

*Cellesse does not reduce the fat in the fatlayer. Therefore women who have severe weight problems or who are obese are discouraged from using Cellesse as the effects are likely to be minimal. This is because the swelling of the fat layer purely by fat is already significant and the effects of a reduction in fluid retention will be minimal in comparison.*

## *When not to use Cellesse*

*Cellesse should not be used by women on areas of the body where they suffer from varicose veins, scars, spots or blemishes, irritated skin or if there is another skin problem or skin disease.*

*In all of these cases, or when in doubt, we suggest you consult your doctor first.*

*In addition, women who have diagnosed heart conditions, diabetes, blood circulation problems or who are pregnant should always consult their doctor before using Cellesse.*

## *Effects during and after use*



*After using Cellesse the skin may become a little red and feel somewhat dry. This is quite normal when massaging and will soon disappear. If needed, some body lotion or cream can be applied after use to soften the skin.*

*Some women may find the use of Celesse slightly painful. As this sensitivity is related to the level of fluidpressure in the fatlayer some women will experience this more than others. During the month this tenderness of the skin can also vary due to the menstrual cycle. In these circumstances we recommend using the lower settings of Celesse untill the skin has become less tender.*

*We also advise women who are prone to bruising of the skin or who have small superficial red veins to be extra careful when using the higher settings, especially in the first weeks of usage. It is also recommended in these cases to use the Vacustop for easier release of the appliance from the skin.*

*In some cases, fine hairs on the legs may be pulled out by the rolling action, if this happens, and causes discomfort, shave the legs before treatment.*

# Instructions for use

## General description

- A** Handgrip
- B** Suction opening (for lifting the skin)
- C** Front roller (motor driven)
- D** Rear roller (motor driven)
- E** Filter (keeps the suction pump free from dust)
- F** On/Off switch and suction control  
*0 = off  
1 = minimum suction power  
5 = maximum suction power*
- G** 'Vacu-Stop' knob (press to switch off the suction pump; then you can easily remove the appliance from your skin)
- H** Transformer plug with power cord (the appliance works on low voltage, so it is safe)
- I** Appliance plug
- J** Connector for appliance plug
- K** Pouch

## Important

- Before connecting the appliance, please check that the voltage indicated on the transformer plug corresponds to the mains voltage in your home.
- Regular cleaning and good maintenance ensures optimum results and guarantees longer life time of Celesse.
- Do not allow the appliance to come into contact with water.
- Store the appliance in its pouch when not in use.
- Keep the appliance away from children.
- Only use your Celesse with the mains cord supplied with the appliance.
- The mains cord of this appliance cannot be replaced. If the cord is damaged, the adaptor must be replaced at the nearest Philips Service Centre.

## *How to get the best from Cellesse*

*Please read the instructions carefully - they are designed to help you get the best possible results.*

*Cellesse has been designed to be used specifically on the upper legs and buttocks where most of the Cellulite is located. We advise not to use Cellesse on the face, the belly or breasts as these areas are more tender.*

## *Frequency of use*

- The ideal frequency of use is every second day. Allow 10-15 minutes for each affected area for the best results.*


## *Intensity*

- Begin each course of treatment on Level 1. During the massage and in subsequent sessions you can gradually increase to higher intensity levels as the skin will get more and more used to the treatment. In this you should establish your own preference as each skin is different.*
- There is no need to apply hard pressure when using Cellesse. Simply let the appliance run over the skin and do the job for you.*


## *Skin Preparation*

- Your skin must be clean, dry and completely free from grease or any other skin preparation.*
- (Cellulite) creams, lotions, talcum powder or any other preparation you normally use can be applied after treatment.*

## *Getting Started*


*Insert the appliance plug into the connector.  
Insert the transformer plug into the wall socket.*







*Set the on/off switch to Level 1 (or the level you are using).*

## *Treatment Techniques*

*Comfort is an important part of Cellesse treatment.  
Try sitting ...*



When you have reached the highest point, press Vacu-Stop and lift the appliance from the skin. The Vacu-Stop is especially useful at the higher intensity levels.



... or standing

Move the appliance slowly and carefully in the right direction, the Vacumotion will lift the skin resulting in the deep intensive massage of the fat layer beneath the skin.

A full movement from just above the knee to the top of the thigh should take around 5 seconds - try counting



Repeat this procedure, treating the entire area with a series of adjacent long, slow, upward strokes.

Gently position the appliance on the skin, a little above the knee, remembering that you are working in a single direction only- towards your heart, i.e. upwards on your legs.

Reposition Celesse next to your original starting point, and begin again.

Find the most comfortable grip and change grip when needed during use.

Establish a set pattern of treatment and follow this routine every time you use Celesse.

After treatment, your skin may feel dry or show a little redness. Don't worry, this is normal when massaging the skin and will soon disappear.


## Cleaning and Maintenance

*Never immerse the appliance in water nor rinse it.*

*For best results, keep Celesse clean at all times, using a damp cloth and, perhaps, a spot of detergent.*


*Do not use alcohol, acetone, gasoline, abrasives, scourers etc. for cleaning the appliance.*

*Pay special attention to the rollers on the underside when cleaning.*

*To clean the rollers, select Level 1, cover the suction opening with your thumb so that the rollers start to rotate, then wipe the rollers clean with a damp cloth, as shown.*

## Storage

*Unplug the appliance after use and store it out of reach of children in the pouch provided.*

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