Always here to help you

Register your product and get support at
www.philips.com/welcome

Question?
Contact
Philips

Sonicare FlexCare Platinum

Quick start guide

Please refer to the user manual for
IMPORTANT SAFEGUARDS

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


PHILIPS
sonicare

YOUR SONICARE (FIG. 1)

- 1 Hygienic travel cap
- 2 Brush head
- 3 Handle
- 4 Power on/off button
- 5 Intensity level indicator
- 6 Intensity level (+ and -)
- 7 Brushing modes
- 8 Mode button
- 9 Battery level indicator
- 10 UV sanitiser and base with integrated charger and cord wrap (specific types only)
- 11 Mains cord with plug
- 12 Deluxe charger (specific types only)
- 13 Charger cover with brush head holders and base with cord wrap (specific types only)
- 14 Travel charger

Note: The contents of the box may vary based on the type purchased.

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
Charge for 24 hours.

Review your user manual for additional instructions.

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
Brushing modes: Before you switch on the Sonicare, press the mode button to toggle between modes.

Intensity levels: To adjust the intensity level to the next higher or lower level, you can press the + button to increase intensity or the - button to decrease intensity at any time during the brushing cycle.

Pressure Sensor: A change in vibration (and a slight change in sound) alerts you when you apply too much pressure while brushing.

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All rights reserved. PHILIPS and the Philips shield are registered trademarks of KPNV. Sonicare and the Sonicare logo are trademarks of Philips Oral Healthcare, Inc. and/or KPNV.