

# PHILIPS

Eye energizing device

ActivEyes

Optimode Technology

Multi-temperature surface

Silicone massage surface

USB cable, Travel pouch



MS3020/10



## Energized, refreshed eyes

One tool to fight puffiness, fatigue, dark circles

98% of consumers agreed that ActivEyes helped their eyes radiate more energy, instantly.\*

### Refresh, relax and stimulate with one device

- Optimode technology addresses 3 main eye concerns
- Refresh: 2-min cool mode reduces puffiness and refreshes
- Relax: 4-min warm mode relieves eye fatigue and relaxes eyes
- Stimulate: 3-step mode combats dark circles under the eyes

### Dual-sided treatment head

- Concave surface for comfortable warm and cool treatments
- Massage surface with soft silicone touchpoints

### Convenient

- Light & sound guides lead you easily through treatments
- Neoprene travel pouch for protection

### Portable

- USB charging for better portability



reddot award 2017  
winner

# Highlights

## Optimode Technology



Optimized settings automatically choose the best temperature to address 3 different eye concerns. A cool treatment reduces puffiness and refreshes. A warm treatment relieves fatigue and relaxes. A combination program with massage combats dark circles.

## Refresh Mode



A 1-minute cool treatment per eye at  $-19^{\circ}\text{C}$  reduces puffiness around the eye and leaves you with a fresher look and feeling.

## Relax Mode



A 2-minute warm treatment per eye at  $-40^{\circ}\text{C}$  relaxes, de-stresses and relieves tension in the muscles. It also triggers the production of essential oils in the eye area, preventing dryness and fatigue caused by extensive computer work.

## Stimulate Mode



A 3-step program increases microcirculation, combats dark circles and gives your eyes a healthy glow. Step 1: A warm treatment for 90 seconds per eye at  $-40^{\circ}\text{C}$  stimulates blood flow, triggering fresh oxygen and nutrient delivery. Step 2: A 30-second manual massage under each eye, with silicone touchpoints gently moving in long strokes from nose to ear, prompts lymph drainage. Step 3: A 15-second cold treatment per eye adds a final refreshing touch. The  $-19^{\circ}\text{C}$  temperature setting constricts blood vessels, making dark circles less visible.

## Multi-temperature surface



Our specially designed multi-temperature surface covers the entire eye, while giving you precisely-timed warm or cool treatments at optimal temperatures. Its rounded, concave shape prevents your sensitive inner eye from receiving too much pressure.

## Massage surface



Our manual massage surface has 60 silicone touchpoints to help you stimulate lymph drainage.

## Light and sound guides



Light and sound guides indicate when the device is ready and lead you through each treatment. A blinking light signals when ActivEyes is ready to use. Blue or red lights indicate cool or warm treatments, respectively. Sound alerts tell you when to switch eyes or move to the next treatment.

## USB cable



ActiveEyes comes with a USB cable, for more flexible charging via your computer – whether you're at work, on the road, or at home.

## Travel pouch

Stylish neoprene pouch protects your device when you travel.



**iF DESIGN AWARD 2017**  
ActivEyesIF DESIGN AWARD 2017



**reddot award 2017**  
winner

**Red Dot Award 2017: Winner**  
ActivEyesRed Dot Award 2017: Winner

# Specifications

## ActivEyes

**Optimode Technology:** Optimal cure per eye concern

**Refresh mode for puffiness:** 2-min cool mode at ~19°C

**Relax mode for eye fatigue:** 4-min warm mode at ~40°C

**Stimulate for dark circles:** 3-step combination program

**Dual-sided treatment head:** Multi-temperature surface, Massage surface

**Treatment head:** Concave for full eye coverage, Precision tip

**Massage surface:** Silicone touchpoints

## Accessories

**Pouch:** Travel pouch

## Ease of use

**3 optimal settings:** Refresh, Relax, Stimulate

**Fits into your daily routine:** Use day or night

**Light guide:** Blinks for preparation time, Blue or red for chosen setting

**Sound guide:** Left eye/right eye treatment, Change treatment surfaces

**Display:** 3 optimal settings, Massage head indicator, Switch eye indicator, Not available indicator, Battery indicator

**Operation:** Cordless use

**Cleaning:** Wipe with damp towel

## Design

**Color:** Black

## Power

**Battery type:** Lithium-ion

**Run time:** Up to 28 treatments\*

**Charging:** USB charging, 2 hours on USB3.0 port\*

**Automatic voltage:** 100-240 V

## Service

**2-year warranty:** On the handle



\* 98% of consumers agreed that ActivEyes helped their eyes radiate more energy, instantly - test conducted with 77 Chinese male respondents by an independent agency in 2016

\* Up to 28 treatments per full charge - number of treatments depend on the selected mode

\* 2 hours full charge on USB3.0 port - charging may take longer if another source is used.