

PHILIPS

AVENT

Nipplette

Single



SCF152/01



A simple solution for inverted nipples

Effective and easy to use

Flat or inverted nipples can be a potential problem for women wishing to breastfeed. Through gentle suction the Nipplette pulls the nipple out into a small plastic thimble-like cup.

Simple solution for inverted nipples

- Clinically proven
- Non-surgical – uses gentle suction
- Simple and discreet

When to use Nipplette

- Use with inverted or flat nipples

Highlights

Non-surgical – uses gentle suction

In a matter of weeks of continued wear, the nipple will stay erect. Ideally use before becoming pregnant. However, you can also use in the first six months of pregnancy to achieve long-lasting correction. Any nipple inversion which is not congenital but has occurred recently should be immediately checked by a

doctor before the Niplette is used. Use is not recommended in the last trimester of pregnancy.

Clinically proven

The Niplette: an instrument for the non-surgical correction of inverted nipples DD McGeorge FRCS (Plast). British Journal of Plastic Surgery, December 1993.

Specifications

Country of origin

England

What is included

Disposable breast pad: 2 pcs

Niplette: 1 pcs

Development stages

Stages: 0 - 6 months, Pregnancy

Design

Discreet design

Ease of use

Easy concealing under clothing: Push syringe to release vacuum

Easy to clean

Clean in hot soapy water

Functions

Correct inverted nipples: Use gentle suction

Material

Breast pads: Dermatologically tested, Natural materials

Silky soft feel

Silky soft topsheet: Breathable natural material

