

Philips
Tubular Fluorescent

20W
T12
Plant and Aquarium

046677392277



Long life

Philips Linear Fluorescent plant and aquarium bulbs have a unique spectral range designed to enhance growth of indoor and aquatic plants. The color of the lamp also enhances the appearance of aquarium fish and decorations.

Light effect

- Enhances the appearance of aquarium fish and decorations
- The right spectrum helps accelerate growth of indoor plants

Long lifespan

- Lasts at least 9000 hours

PHILIPS

Tubular Fluorescent
20W T12, Plant and Aquarium

Specifications

Technical specifications

- Lamp holder: T12
- Placement: Aquarium, Fluorescent fixture

- Indoor
- Number of lamps: 1
- Wattage: 20 W

046677392277

Highlights

For freshwater aquariums

Enhances the appearance of aquarium fish and decorations

For growing indoor plants

The right spectrum helps accelerate growth of indoor plants

Lasts at least 9000 hours

Lasts at least 9000 hours



Issue date 2016-01-29

Version: 5.0.0

12 NC: 9278 830 00700
UPC: 0 46677 39227 7

© 2016 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Average life under specified test conditions with lamps turned off and restarted no more frequently than once every 3 operating hours. Lamp life is appreciably longer if lamps are started less frequently.