



Highlight your elegance

Add a touch of luxury and brilliance to your home with Philips myLiving Treasure LED Suspension light in gold glass shade and premium faceted glass crystals. Experience warm white or cool daylight effects at a simple flick of the switch.

Match your interior style

- Coordinated range

Premium quality

- Premium materials

Advanced

- Built-in VariLume technology for easy brightness control

Light effect

- Multiple light sources

Philips myLiving
Suspension light

36899 Treasure

gold
LED

myLiving



36899/10/C5

PHILIPS

Suspension light
36899 Treasure gold, LED

Specifications

Design and finishing

- Material: glass
- Color: gold

Product dimensions & weight

- Height minimum: 35 cm
- Height maximum: 150 cm
- Length: 40 cm
- Width: 40 cm
- Net weight: 4.900 kg

Technical specifications

- Mains power: Range 220 V - 240 V
- Number of bulbs: 1
- Wattage bulb included: 36 W
- Total lumen output fixture: 1500 lm
- Fixture dimmable: No
- LED

- IP code: IP20, protection against objects bigger than 12.5 mm, no protection against water
- Class of protection: I - earthed

Service

- Warranty: 1 year(s)

Packaging dimensions & weight

- Height: 45 cm
- Length: 47.5 cm
- Width: 47.5 cm
- Weight: 7.800 kg

Miscellaneous

- Especially designed for: Living- & Bedroom
- Style: Historical
- Type: Suspension light

3689910C5

Highlights

Coordinated range



Complete range of coordinated designs available for you to choose.

Premium materials



Made of premium materials to assure the quality of the product.

Multiple light source

Multiple light sources

VariLume technology



Advanced electronic controller, for variable light levels - at the flick of a normal switch



Issue date 2018-07-03

Version: 1.1.1

12 NC: 9150 052 82801
EAN: 69 47830 43099 4

© 2018 Philips Lighting Holding B.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Philips Lighting Holding B.V. or their respective owners.

www.philips.com