



Philips myLiving
Ceiling light

50166 Krilyst

white
LED

myLiving



50166/11/66

Highlight your elegance

Bask in the glorious sparkle of the Philips myLiving Krilyst luminaire. Enjoy the look of its crystals hanging elegantly from the ceiling and control the brightness of the LED light with the wall switch. Perfect for livingrooms or bedrooms.

Premium quality

- Premium materials

Energy saving

- Built-in LED, as part of the system

PHILIPS

Ceiling light
50166 Krilyst white, LED

Specifications

Design and finishing

- Color: white

Product dimensions & weight

- Height: 23 cm
- Length: 55 cm
- Width: 55 cm
- Net weight: 4.200 kg

Technical specifications

- Mains power: Range 100 V - 240 V, 50-60 Hz
- Number of bulbs: 12
- Wattage bulb included: 0.5 W
- Maximum wattage replacement bulb: 30 W
- Fixture dimmable: No
- LED

- Built in LED
- IP code: IP20, protection against objects bigger than 12.5 mm, no protection against water
- Class of protection: I - earthed

Packaging dimensions & weight

- Height: 21 cm
- Length: 71.5 cm
- Width: 71.5 cm
- Weight: 6.800 kg

Miscellaneous

- Especially designed for: Living- & Bedroom
- Style: Historical
- Type: Ceiling light

501661166

Highlights

Premium materials



Made of premium materials to assure the quality of the product.

Integrated LED light source



A light source you can trust. With integrated LED light source, you won't have to worry about maintenance or lamp replacement while having the perfect light ambiance in your home.



Issue date 2018-03-16

Version: 1.0.1

12 NC: 9150 051 37001
EAN: 69 47830 42779 6

© 2018 Philips Lighting Holding B.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Philips Lighting Holding B.V. or their respective owners.

www.philips.com