

# PHILIPS

Cordless Wet & Dry epilator

Satinelle Prestige

For legs, body and face

4 in 1 Epilator + Shaver with Body Care

Extra wide epilation head

Rechargeable + 40 min of use



BRE648/00



## More hair removal power + tweezing actions in a Philips epilator

with an award winning design and Body Care routines

Satinelle Prestige 4 in 1 Epilator & Shaver with Body Care is a complete hair removal kit designed to effectively remove hair on the whole body - even sensitive areas - for weeks of satin smooth skin. Special ceramic tweezing disks rotate quickly and capture more thick and thin hairs.

### Exceptional results

- Special ceramic tweezing discs grip more hairs
- Extra wide epilator head

### Easy and effortless usage

- Award-winning design\*
- Cordless wet and dry for use in bath or shower
- Opti-light helps you target and remove even trickiest hairs
- First epilator with S-shaped handle

### Personalized body and face care routines

- Body exfoliation brush removes dead skin cells
- Body massager to relax and for a radiant looking skin
- Shaving head and trimming comb for a close shave
- Facial area and delicate area cap to remove unwanted hairs
- Includes massage cap
- Includes skin stretcher


# Highlights

## Our fastest epilation ever



Our epilator head is unique in being made from a textured ceramic surface that gently extract even the finest hairs and those shorter than wax. Now with faster disc rotation than ever before (2200RPM) for our fastest hair removal.

## Extra wide epilator head



Extra wide epilation head grabs more hairs on the first pass for faster hair removal without pulling sensitive skin.

## Award-winning design\*



Award-winning design\* for effortless hair removal

## Cordless wet and dry usage



Designed with an anti-slip grip, ideal for use with water. Enables a more comfortable, gentle experience in your shower or bath. You can use it cordless for best convenience.

## Opti-Light



Opti-light helps you target and remove even trickiest hairs

## Body exfoliation brush



Our body exfoliation brush removes dead skin cells and helps to prevent ingrown hairs. The hypo-allergenic 48.200 fine bristles gently yet effectively sweep away dead skin cells & stimulate the regeneration of the skin surface. This device will exfoliate your skin more effectively than manual treatment alone.

## Body massager



Our body massager gives you a relaxing treatment so you can enjoy radiant-looking skin.

## Shaving head and trimming comb



The shaving head provides a close shave and more gentleness in different body areas. It comes with a trimming comb to style your bikini area.

## Facial and delicate area cap



For more gentleness in all different body areas it includes a facial area cap to easily remove unwanted facial hairs and a delicate area cap for hair from underarm and bikini

# Specifications

## Accessories

Shaving head  
Trimming comb  
Body Exfoliation Brush  
Body Massager  
Delicate area cap  
Skin stretcher cap  
Facial area cap  
Massage cap  
Pouch: Basic pouch  
Cleaning brush

## Ease of use

Wet and dry use  
Opti-light  
Cordless  
Handle: S-shape handle

## Performance

Epilation system: Patented epilation system  
Epilation discs: Ceramic discs

## Features

Speed settings: 2 settings

## Technical specifications

Number of catching points: 32  
Number of discs: 17  
Voltage: 15V / 5.4W  
Tweezing action speed 1: 64000 per minute  
Tweezing action speed 2: 70400 per minute

## Power

Quick charge  
Battery Type: Lithium-ion  
Charging: Rechargeable, 1.5 hour charging time  
Usage time: up to 40 minutes

