

# PHILIPS

## Milk circuit cleaner sachets

Same as CA6705/60

For 6 uses - Use monthly

Prolong machine lifetime

Brew better tasting coffee



CA6705/10

## Brew better tasting cappuccino

### Keep milk circuit hygienic and clean

With the Philips milk circuit cleaner you can sanitize the milk circuits of your espresso machine or milk frother. It ensures effective removal of milk residue from all components. We recommend to carry out this cycle once per month.

#### **Sensational milk foam quality**

- Improve the taste of your coffee specialty

#### **Effective protection of your espresso appliance**

- Protects the milk circuits from milk residue clogging

#### **Prolong lifetime**

- Regular cleaning prolongs lifespan of milk frothing parts

#### **Reliable machine protection**

- Only use Philips consumables

# Highlights

## Perfect coffee taste




Thanks to the milk circuit cleaner cleaning solution you can effectively remove milk residues that might affect milk taste.

## Protects the milk circuits



Protects the milk circuits from milk residue clogging

## Prolongs machine lifetime



Regular cleaning prolongs the life of all milk frothing parts.

## Philips Original



Only use Philips consumables to make sure your machine runs smoothly for longer. Philips consumables are the only recommended consumables for Philips and Saeco machines.

# Specifications

## Country of origin

Made in: Germany

## Technical specifications

Includes: 6x milk circuit cleaner sachets

## Weight and dimensions

Weight of product: 0.1 kg

Quantity: 6 tablets of 1.6 g

