

PHILIPS

Quick Clean Pod Cartridge

2-pack

Up to 6 mos hygienic shaving

Quick Clean Pod compatible



CC12/51



Enjoy hygienic shaving every day 10 times more effective than cleaning with water*

Using the Quick Clean Pod after every shave will keep your shaver like new, removing cut hair effectively. It leaves your shaver 10 times cleaner versus using water alone*.

Deep cleaning in just 1 minute

- Philips exclusive hygienic formula + multi-grid filter
- A cartridge is effective for up to three months
- Ensures a hygienic shave for up to 6 months
- Quick Clean Pod compatibility

Active lubrication

- Lubricants keep your shaver performing at its best

Refreshing scent

- Fresh fragrance for a clean-feeling shave

Skin-friendly formula

- Completely alcohol-free

Highlights

Thorough cleaning



Philips exclusive hygienic formula combined with a multi-grid filter eliminates cut hair effectively and leaves your shaver 10 times cleaner versus using water alone* and prevent 99.9%** of bacteria growth.

Active lubrication



The formula, enriched with active lubricants, protects your shaving heads from friction and wear, keeping your shaver performing at its best for longer.

Refreshing scent



The cleaning cartridge's special fragrance smells fresh and gives you a clean feel while shaving.

Skin-friendly formula



Completely alcohol-free, the cleaning fluid is skin-friendly and designed to offer protection from skin irritation.

Long-lasting



Each Philips cleaning cartridge is effective for about 30 cleaning cycles with daily usage and up to three months with weekly usage. That's an entire season's-worth of clean, hygienic shaving.

6-month supply



Enjoy up to 6 months of fresh shaving with this cartridge 2-pack.

For Quick Clean Pod



The cartridge is compatible with the Philips Quick Clean Pod.

Specifications

Items included

Cartridge: 2 pcs

Capacity

Cartridge: 2x 5.41 fl oz / 2x 160 ml

