

# PHILIPS

## Quick Clean Pod Cartridge

5+1-pack

Up to 18 mos hygienic shaving

Quick Clean Pod compatible



CC16/51



## Enjoy hygienic shaving every day 10 times more effective than cleaning with water\*

Using the Quick Clean Pod after every shave will keep your shaver like new, removing cut hair effectively. It leaves your shaver 10 times cleaner versus using water alone\*.

### Deep cleaning in just 1 minute

- Philips exclusive hygienic formula + multi-grid filter
- A cartridge is effective for up to three months
- Ensures a hygienic shave for up to 18 months
- Quick Clean Pod compatibility

### Active lubrication

- Lubricants keep your shaver performing at its best

### Refreshing scent

- Fresh fragrance for a clean-feeling shave

### Skin-friendly formula

- Completely alcohol-free

# Highlights

## Thorough cleaning



Philips exclusive hygienic formula combined with a multi-grid filter eliminates cut hair effectively and leaves your shaver 10 times cleaner versus using water alone\*

## Active lubrication



The formula, enriched with active lubricants, protects your shaving heads from friction and wear, keeping your shaver performing at its best for longer.

## Refreshing scent



The cleaning cartridge's special fragrance smells fresh and gives you a clean feel while shaving.

## Skin-friendly formula



Completely alcohol-free, the cleaning fluid is skin-friendly and designed to offer protection from skin irritation.

## Long-lasting



Each Philips cleaning cartridge is effective for about 30 cleaning cycles with daily usage and up to three months with weekly usage. That's an entire season's-worth of clean, hygienic shaving.

## 18-month supply



Enjoy up to 18 months of fresh shaving with this cartridge 6-pack.

## For Quick Clean Pod



The cartridge is compatible with the Philips Quick Clean Pod.

# Specifications

**Items included**

Cartridge: 6 pcs

**Capacity**

Cartridge: 6x 5.41 fl oz / 6x 160 ml

