

PHILIPS

Lid

Stainless steel

CP0348



to replace your current removable lid

this is the removable lid for the QC Basket

Ease of using Philips spare parts

- Easy renew your product with original Philips parts

Highlights

Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

Specifications

Replaceable part

Fits product type:: HD9620/0x, HD9621/1x, HD9621/2x, HD9621/9x, HD9622/2x, HD9622/3x, HD9622/4x, HD9622/9x, HD9623/1x, HD9627/1x, HD9627/2x, HD9621/4x

