

# PHILIPS

Juice collector

black

CP0439



## to replace your current juice collector

Philips' unique pre-clean function allows you to quickly rinse clean your juicer between different juices or after juicing. By pouring water into the pusher, you can create a water fountain, rinsing away unwanted leftovers.

### **Ease of using Philips spare parts**

- Easily renew your product with original Philips parts

# Highlights

## Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

# Specifications

## Replaceable part

**Fits product type::** HR1916/xx, HR1919/xx, HR1920/xx, HR1921/xx, HR1922/xx, HR1925/xx

