

PHILIPS

Complete milk carafe

Exprelia Evo

Quick clean after use

Thorough clean weekly

Not dishwasher proof

CP0502



Complete carafe kit

Please check specifications for compatible product

This complete milk carafe, including milk tube, makes your milk deliciously foamy. The carafe can also be stored conveniently in the fridge. Compatible with Exprelia Evo.

Ease of using Philips spare parts

- Easy renew your product with original Philips parts

Highlights

Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

Specifications

Suitable for:

Exprelia series: HD8855/01, HD8855/09, HD8857/01, HD8857/09

