

# PHILIPS

Holder for shaving  
device

1 piece

CP9301



## A brace for fastening your shaving heads

Go to the Specs Tab to check your matching product

The shaving head holder keeps the heads of your shaving device in place. Is your current one lost or broken and the rest of your shaving device still perfect? Order a new holder here!

### **Ease of using Philips spare parts**

- Easily renew your product with original Philips parts

# Highlights

## Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

# Specifications

## Replaceable part

**Fits product types:** PT710, PT725, PT735, PT736, PT860, PT870, PT875, PT876

**Fits product types::** PT710/14, PT710/15, PT710/16, PT710/17, PT710/20, PT710/41, PT715/14, PT715/15, PT715/17, PT715/41, PT720/14, PT720/15, PT720/16, PT720/17, PT720/41, PT725/15, PT725/17, PT725/20, PT730/14, PT730/15, PT730/17, PT730/20, PT730/41, PT735/15, PT735/16, PT735/17, PT735/20, PT735/41, PT736/18

**Fits product types::** PT710, PT715, PT720, PT725, PT730, PT735, PT736

**Fits product type:** PT711, PT717, PT721, PT722, PT723, PT726, PT728, PT729, PT731, PT732, PT734, PT739, PT761, PT762, PT786, PT845, PT861, PT862, PT868, PT871, PT872, PT874, PT878, S1100, S1110, S1150, S1310, S1510, S1520, S1560, S1570, S3110, S3120, S3510, S3520

**Fits product types::** PT711/16, PT711/20, PT717/16, PT717/17, PT720/49, PT721/16, PT723/16, PT723/17, PT723/20, PT723/80, PT723/81, PT727/16, PT727/17, PT727/20, PT727/80, PT731/16, PT731/17, PT737/16, PT737/17, PT737/18, PT739/14, PT739/18, PT739/20

**Fits product type::** PT720, PT856, PT860, PT867

