

PHILIPS

Detachable water
tank for your iron

Wardrobe Care

Blue

CRP164



Water tank for your Wardrobe Care

Check compatibility data

Want to extend the lifetime of your Philips Wardrobe Care? The cartridge water tank is designed for use with anti-calc cartridges to prevent limescale from building up inside the system.

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types: GC9920/07, GC9920/27

Fits product types:: To be defined

Fits product type:: GC9940/05, GC9940/07, GC9940/27, GC9950/07

