

# PHILIPS

Detachable water  
tank for your iron

PerfectCare

Transparent purple

CRP174



## Extra large water tank for your Wardrobe Care

Check compatibility data

This extra-large, 1.5-litre water tank allows up to two hours of ironing with your PerfectCare iron, without refilling. Easily replaceable!

### **Ease of using Philips spare parts**

- Easily renew your product with original Philips parts

## Highlights

### Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

## Specifications

### Replaceable part

**Fits product types::** To be defined

**Fits product type::** GC9240/02, GC9240/07, GC9240/27

