



Detachable water tank for your iron

- PerfectCare
- Transparent purple

CRP174

Extra large water tank for your Wardrobe Care

Check compatibility data

This extra-large, 1.5-litre water tank allows up to two hours of ironing with your PerfectCare iron, without refilling. Easily replaceable!

Benefits

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Features

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Replaceable part

Fits product types: To be defined

Fits product type: GC9240/02
GC9240/07

PHILIPS

GC9240/27



data subject to change
2019, September 10

Version: 2.0.1
EAN: 08710103564928

© 2019 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com