

PHILIPS

Basket

3 in 1

white

220 V

CRP231



Basic part for steriliser

Find matching products on the specifications tab

The top basket of your electric steam steriliser is meant to sterilise smaller baby feeding items, like soothers, bottle teats, bottle screw rings and bottle dome caps.

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types:: SCF284/02, SCF284/03, SCF285/01, SCF285/02

