



Basket

- 3 in 1
- white
- 220 V

CRP231

Basic part for steriliser

Find matching products on the specifications tab

The top basket of your electric steam steriliser is meant to sterilise smaller baby feeding items, like soothers, bottle teats, bottle screw rings and bottle dome caps.

Benefits

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Features

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Replaceable part

Fits product types: SCF284/02

SCF284/03
SCF285/01

PHILIPS

SCF285/02



data subject to change
2019, April 10
Version: 4.0.1
EAN: 08710103564331

© 2019 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com