

PHILIPS

Basket

4 in 1

White/grey

220 V

CRP235



Basic part for steriliser

Find matching products on the specifications tab

This basket for your electric steam sterilizer perfectly fits all baby feeding items you want to sterilize: both wide and narrow neck Avent bottles, parts of your breast pump and other accessories.

Ease of using Philips spare parts

- Easy renew your product with original Philips parts

Highlights

Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types:: SCF286/01, SCF286/02, SCF286/03, SCF286/05, SCF287/01, SCF287/02

