

PHILIPS

Trimmer

For QG3190 and 3193

CRP295/01



Shaving unit with foil for body groomer

Check compatibility data

This trimmer, part of your body groom set, is designed to catch all types of body hair with one-stroke efficiency. Guarantees an excellent and fast result!

Ease of using Philips spare parts

- Easy renew your product with original Philips parts

Highlights

Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types:: QG3190, QG3193

