

PHILIPS

Trimmer

For QG3190 and 3193

CRP295/01



Shaving unit with foil for body groomer

Check compatibility data

This trimmer, which is part of your body groom set, is designed to efficiently catch all types of body hair with one-stroke. Guarantees excellent and fast results!

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types:: QG3190, QG3193

