



Trimmer

• For QG3190 and 3193

CRP295/01

Shaving unit with foil for body groomer

Check compatibility data

This trimmer, which is part of your body groom set, is designed to efficiently catch all types of body hair with one-stroke. Guarantees excellent and fast results!

Benefits

[Ease of using Philips spare parts](#)

- Easily renew your product with original Philips parts

Features

[Easy replacement of product parts](#)

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

PHILIPS

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Replaceable part

QG3193

Fits product types: QG3190



data subject to change
2018, December 31

Version: 2.0.1
EAN: 08710103568599

© 2018 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com