


PHILIPS

Beard trimmer comb

For QG3190 and 3193

Small comb

CRP296/01



Accessory for your grooming kit

Check compatibility data

With this comb, which is part of your grooming kit, you can cut and style your beard just the way you like it. Guarantees unprecedented precision to create exactly the look you had in mind!

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types:: QG3190, QG3193

