

# PHILIPS

Shaving head holder

For HQ8445-8899

CRP343/01



## A brace for fastening your shaving heads

Check compatibility data

The shaving head holder keeps the heads of your shaving device in place. Is your current shaving head lost or broken, but the rest of your shaving device still works perfectly? Order a new holder here!

### **Ease of using Philips spare parts**

- Easily renew your product with original Philips parts

# Highlights

## Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

# Specifications

## Replaceable part

**Fits product types::** HQ8445, HQ8825, HQ8830, HQ8845, HQ8850, HQ8865, HQ8870, HQ8875, HQ8880, HQ8882, HQ8885, HQ8890, HQ8893, HQ8894

**Fits product types:** 8831XL, 8845XL, 8846XL, 8867XL, 8881XL, 8892XL, 8894XL

