

**PHILIPS**

**AVENT**

Lid and clips for  
Microwave Steriliser

microwave steriliser

CRP408



## To cover your Microwave Steriliser

Find matching products on the specifications tab

The Philips Avent Microwave Steam Sterilizer is for sterilizing bottles and other means you use for baby feeding. This lid is meant to cover it. Keep it unopened after sterilizing and the sterilizer's content remains sterile for 24 hours.

### **Ease of using Philips spare parts**

- Easy renew your product with original Philips parts

# Highlights

## Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

# Specifications

## Replaceable part

**Fits product types ::** SCD279/00, SCD283/00, SCF271/01, SCF271/02, SCF271/06, SCF271/07, SCF271/20, SCF271/42, SCF271/51, SCF277/01, SCF277/51

