

PHILIPS

AVENT

Lid and clips for
Microwave Steriliser

microwave steriliser

CRP408



To cover your Microwave Steriliser

Find matching products on the specifications tab

The Philips Avent Microwave Steam Steriliser is for sterilising bottles and other accessories you use for baby feeding. This lid is used to cover it. Keep it unopened after sterilising and the steriliser's content remains sterile for 24 hours.

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types:: SCD279/00, SCD283/00, SCF271/01, SCF271/02, SCF271/06, SCF271/07, SCF271/20, SCF271/42, SCF271/51, SCF277/01, SCF277/51

