

PHILIPS

AVENT

Ice packs

out and about

CRP409



To keep expressed milk cool

Find matching products on the specifications tab

This flexible gel pack is an essential part of the Avent Deluxe Bag. It allows you to keep your breast milk cool for 4 hours outside the refrigerator. (Note: 1 pack per order)

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types:: SCF290/11, SCF290/12, SCF290/13, SCF290/15, SCF290/16, SCF290/20, SCF290/28, SCF290/60, SCF292/01, SCF292/13, SCF292/15, SCF292/31, SCF294/00, SCF294/02, SCF294/21, SCF300/08, SCF300/12, SCF300/13, SCF300/15, SCF300/20, SCF300/60, SCF302/01, SCF302/13, SCF302/60, SCF304/02, SCF304/60, SCF310/12, SCF310/13, SCF310/20, SCF310/60, SCF312/01, SCF312/13, SCF312/60, SCF314/02

Fits product type:: SCF290/98, SCF290/99

