

PHILIPS

Detachable water
tank for your iron

GC9000 series

Blue

CRP589



For storing the water in your iron

Check compatibility data

Filling this detachable water tank with water is necessary to generate steam with your iron. The container allows you to iron for two hours without needing to refill.

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types:: GC9020/02, GC9020/07, GC9040/02, GC9140/02, GC9140/07

