

PHILIPS

Cappuccinatore

for Milkfrother

CRP993



For the automatic frothing of milk

Please check specifications for compatible products

Do you like a good cappuccino or latte macchiato? This black cappuccinatore or milk frother makes your milk deliciously foamy. It fits different espresso machines.

Ease of using Philips spare parts

- Easy renew your product with original Philips parts

Highlights

Easy renew of product parts

From time to time your product needs a facelift but it was never so easy like with Philips consumer replacement parts to renew your product. This all with guaranteed Philips quality.

Specifications

Replaceable part

Fits product types: 10003804, 10003805

Fits product types :: 10003973

