

PHILIPS

ZOOM!

Take-home
whitening treatment

NiteWhite

6 syringes

18 applications

in several concentrations

DIS114/00



Whiter nights

Professional take-home whitening treatment

With Philips Zoom NiteWhite, your most confident, radiant smile isn't just whiter, but healthier too. And because you whiten at home in the evening, a whiter smile has never been more convenient.

A whiter, healthier smile

- Healthier, stronger enamel
- Your smile can be noticeably whiter

Professional-quality results

- 14 million smiles and counting****
- Proven safe and effective on teeth and gums
- Recommended and overseen by your dental professional

Whitening on your time

- The perfect fit for your smile
- Whiten at home on your schedule

Highlights


Whiten at home

Zoom NiteWhite lets you whiten your teeth when it's most convenient for you. Just wear your trays for either one to two hours a day or overnight, depending on the formula your dentist recommends. And don't worry if you miss a day. You can start and stop your whitening program to accommodate your busy schedule.

Healthier enamel

Philips Zoom NiteWhite is a take-home whitener with the power of three: Creating a lustrous smile, smoothing and shining teeth, and managing sensitivity. That means it's making your smile healthier as it whitens by improving your enamel, which will be smoother and more lustrous by the time you finish treatment.

Noticeably whiter teeth



With Philips Zoom NiteWhite, your whiter smile can be noticeably whiter fast, with maximum results in one to two weeks.

The perfect fit

With Philips Zoom NiteWhite, your dental professional makes a custom whitening tray to perfectly fit the contours of your teeth. That means you get the same bright results across your entire smile!

Proven safe and effective

Philips Zoom NiteWhite is backed by studies that have shown it is safe and effective for your teeth and gums. And because your dental professional recommends it, you can be sure that your results will be consistent and safe on any fillings or restorations you might have.

Recommended professionally

Philips Zoom NiteWhite is recommended and overseen by your dental professional, so you know it's the best, healthiest option for your teeth.

Chosen by millions

Philips Zoom NiteWhite is part of the popular Philips Zoom whitening family, the #1 patient-requested professional whitening brand that's been trusted by over 14 million people already****. It also comes from the same creators of popular brands you know and love, including Philips Sonicare and BreathRx, and was designed to help improve the health of your smile.

Specifications

Performance

Whitening results: Maximum results in 7-14 days, Whiter teeth in 3 days

Proven safe: With little to no sensitivity

Dual-barrel syringe: No refrigeration required

Customizable

Ingredients

Carbamide Peroxide: 22% (equiv. to 7.96% HP)

ACP: Amorphous calcium phosphate

Flavor: Mint

Items included

Tray case: Holds custom whitening tray

Spa Bag: Yes*

