

# PHILIPS

## ZOOM!

Take-home  
whitening treatment

DayWhite

3 syringes

9 applications

9.5% Hydrogen Peroxide

DIS120/00



## Brighter days

### Professional take-home whitening treatment

With Philips Zoom DayWhite, your most confident, radiant smile isn't just whiter, but healthier too. And because you whiten on your time at home, that brighter smile has never been more convenient.

#### **A whiter, healthier smile**

- Healthier, stronger enamel
- Your smile can be noticeably whiter

#### **Professional-quality results**

- 14 millions smiles and counting\*\*\*\*
- Proven safe and effective on teeth and gums
- Recommended and overseen by your dental professional

#### **Whitening on your time**


- The perfect fit for your smile
- Whiten at home on your schedule

# Highlights

## Whiten at home

Philips Zoom DayWhite lets you whiten your teeth when it's most convenient for you. Just wear your trays for 15 to 30 minutes twice daily, depending on the formula your dentist recommends. And don't worry if you miss a day. You can start and stop your whitening program to accommodate your busy schedule.

## Noticeably whiter teeth



With Philips Zoom DayWhite, your whiter smile can be noticeably whiter fast, with maximum results in just seven to 14 days.

## Healthier enamel

Philips Zoom DayWhite is a take-home whitener with the power of three: Creating a lustrous smile, smoothing and shining teeth, and managing sensitivity. That means it's making your smile healthier as it whitens by improving your enamel, which will be smoother and more lustrous by the time you finish treatment.

## The perfect fit

With Philips Zoom DayWhite, your dental professional makes a custom whitening tray to perfectly fit the contours of your teeth. That means you get the same bright results across your entire smile!

## Proven safe and effective

Philips Zoom DayWhite is backed by studies that have shown it is safe and effective for your teeth and gums. And because your dental professional recommends it, you can be sure that your results will be consistent and safe on any fillings or restorations you might have.

## Recommended professionally

Philips Zoom DayWhite is recommended and overseen by your dental professional, so you know it's the best, healthiest option for your teeth.

## Chosen by millions

Philips Zoom DayWhite is part of the popular Philips Zoom whitening family, the #1 patient-requested professional whitening brand that's been trusted by over 14 million people already\*\*\*\*. It also comes from the same creators of popular brands you know and love, including Philips Sonicare and BreathRx, and was designed to help improve the health of your smile.

# Specifications

## Performance

**Whitening results:** Maximum results in 7-14 days, Whiter teeth in 3 days

**Proven safe:** With little to no sensitivity

**Dual-barrel syringe:** No refrigeration required

**Customizable**

## Power of three

**ACP:** Creates a lustrous smile, Smooth and shiny teeth

**Potassium nitrate:** Manages sensitivity

**Fluoride:** For a healthy white smile

## Ingredients

**Hydrogen Peroxide:** 9.5%

**Flavor:** Mint

**ACP:** Amorphous calcium phosphate

## Items included

**Dual-barrel syringes:** 3

**Tray case:** Holds custom whitening tray

**Spa Bag:** Yes\*

