

PHILIPS

Ear thermometer

Accurate temperature reading

Fast measurement in 2 seconds

Syncs with Philips health app



DL8740/15



Connected temperature tracking

Syncs with free Philips HealthSuite health app

The connected ear thermometer allows you to quickly and accurately take readings and store them directly on your smart device using the free Philips HealthSuite health app, so you can conveniently monitor your temperature over time.

Delivers fast and accurate results

- Provides accurate temperature measurements
- Accurate temperature reading in 2 seconds

Personalised tracking in the Philips health app

- Automatically syncs to the Philips health app via Bluetooth®
- Monitor your temperature over time

Comfortable and easy to use

- Suitable for people of all ages
- Gives results in Celsius and Fahrenheit
- Red LED display ring for rapid high-temperature indication



asimpleswitch.com

Highlights

Accurate




Fast and accurate in-ear measurements. The Philips ear thermometer uses in-ear measurement to offer precise readings you can be confident are accurate. The thermometer offers accurate measurements of $\pm 0.2^{\circ}\text{C}$ ($\pm 0.4^{\circ}\text{F}$) within the range of 32.4°C - 42.9°C (90.3°F - 109.3°F).

Phone compatibility



Automatically transfer your measurements to the Philips HealthSuite health app on your smartphone so you can monitor changes over time. Compatible with iOS 9.3 or higher/Android 4.4 KitKat or higher. For an overview of (in)compatible devices, please go to www.philips.com/personalhealthcompatibility

Fast



Specially designed for speed and ease. Just one press of the button and you will have a reading in 2 seconds.

Easy and comfortable to use




The thermometer has a simple one-touch button and a probe tip that is comfortable to insert and suitable for use by all ages (newborn and up). No replacement caps required. Just clean the thermometer tip before use.

Celsius and Fahrenheit



No matter whether you measure in Celsius or Fahrenheit, the Philips HealthSuite Ear thermometer will give you the results you need in the scale you understand

Track your results over time



Determine your normal body core temperature so you can track and spot changes quickly and easily in the Philips HealthSuite health app.

Red LED display



The Philips ear thermometer has a clear LCD display so you can instantly see if the temperature reading is at a normal* or raised level. The light ring turns red for temperatures of 38.0°C (100.4°F) and higher. *Normal temperature ranges differ by age group. 38.0°C (100.4°F) is still within the normal temperature range for babies.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Syncing

Phone compatibility: iOS 8 or higher / Android 4.4 or higher

Data transfer technology: Bluetooth®

Design

Colour: White

Dimensions: 136(L) x 34.5(W) x 57(D) mm

Net. Weight: About 54.6 g

Ease of use

Fast: 2-second measurement

Switch to Celsius and Fahrenheit

Clear LCD display

Probe tip for all ages

Can be used as stand-alone

Technical

Accuracy: ± 0.2 °C (± 0.4 °F)

Measurement range: 32.4 °C–42.9 °C (90.3 °F–109.3 °F)

Power

Power source: 6 V (2 x AAA alkaline batteries)

Accessories

2 x AAA batteries

Service

2 year guarantee

