

# PHILIPS

## Upper arm blood pressure monitor

Blood pressure measurement

Heart rate measurement

Syncs with Philips health app



DL8760/01



## Connected blood pressure tracking

### Syncs with free Philips HealthSuite health app

The connected upper arm blood pressure monitor allows you to measure your blood pressure & heart rate. Syncs automatically to the free Philips HealthSuite health app so you can conveniently track your measurements over time.

#### **Comfortable and easy to use**

- Large backlit screen clearly shows your measurement results
- No batteries required. Rechargeable using DC charger
- Allows 2 different users to monitor & track their readings

#### **Measure & monitor**


- Measures systolic & diastolic pressure and heart rate
- Classifies your measurement according to the WHO system

#### **Personalized tracking in the Philips health app**

- Automatically syncs to the Philips health app via Bluetooth®
- Monitor your blood pressure over time

# Highlights

## Measure & monitor




Measures systolic & diastolic pressure and heart rate

## Blood pressure classification



Instantly compare your blood pressure against the internationally recognized classifications by the World Health Organization (WHO).

## Track your results over time



Shows your blood pressure and heart rate measurements through an easy-to-use dashboard on the Philips HealthSuite health app, so you can monitor changes over time.

## Phone compatibility




Automatically transfer your measurements to the Philips HealthSuite health app on your smart phone, so you can monitor changes over time. Compatible with iOS 9.3 or higher/Android 4.4 KitKat or higher. For an overview of (in)compatible devices please have a look at [www.philips.com/personalhealthcompatibility](http://www.philips.com/personalhealthcompatibility)

## Supports 2 users




Allows 2 different users to measure their blood pressure levels and heart rate & monitor changes over time in the Philips health app.

## Rechargeable



No batteries required. Rechargeable using DC charger

## Easy-to-read display



Large backlit screen clearly shows your measurement results.

# Specifications

## Design

Color: Black

Device dimensions: 135x 76x 30 mm

Display: LCD display, white backlight

Display dimensions: 87x24 mm

Net. Weight: 265 g

## Ease of use

Can be used as stand-alone

Auto sync data

Blood pressure classification

Irregular heartbeat detector

Measure during inflation

LCD display

Supports 2 users

## Accessories

Cuff

DC charger

Arm blood pressure monitor

## Power

Adaptor type: 6V 1A DC charger

Battery life: 20 days\*

Battery Type: Rechargeable li-polymer batt.

Charging time:  $\pm$  2 hr

## Service

2-year guarantee

## Syncing

Data transfer technology: Bluetooth®

Phone compatibility: iOS9.3+

(iPhone5+)/Android4.4+

## Technical

Measurement pressure: 5.3kPa-30.7kPa (40-230mmHg)

Heart rate accuracy:  $\pm$ 5% of result on display

Pulse value measurement range: 40-199

heart beats/minute

Rated cuff pressure: 0kPa-40kPa (0mmHg-300mmHg)

Cuff size of the upper arm: 22 ~42 cm

Pressure accuracy: 5-40°C within  $\pm$ 0.4kPa (3mmHg)

Operating environment: Atmospheric pressure:86-106KPa, Relative humidity: < 85%RH, Temperature: 5°C to 40°C

