

# PHILIPS

## Wrist blood pressure monitor

Blood pressure measurement

Heart rate measurement

Syncs with Philips health app



DL8765/15



## Connected blood pressure tracking

### Syncs with free Philips HealthSuite health app

The connected wrist blood pressure monitor allows you to measure your blood pressure and heart rate. Syncs automatically to the free Philips HealthSuite health app so you can conveniently track your measurements over time.

#### Comfortable and easy to use

- Portable wrist unit to measure your blood pressure anywhere
- Large backlit screen clearly shows your measurement results
- No batteries required: rechargeable using USB cable
- Allows two different users to monitor and track their readings

#### Measure and monitor


- Measures systolic and diastolic pressure and heart rate
- Classifies your measurement according to the WHO system

#### Personalised tracking in the Philips health app

- Automatically syncs to the Philips health app via Bluetooth®
- Monitor your blood pressure over time

# Highlights

## Measure and monitor




Measures systolic and diastolic pressure and heart rate

## Blood pressure classification



Instantly compare your blood pressure against the internationally recognised classifications by the World Health Organisation (WHO).

## Track your results over time



Shows your blood pressure and heart rate measurements through an easy-to-use dashboard on the Philips HealthSuite health app, so you can monitor changes over time.

## Phone compatibility



Automatically transfer your measurements to the Philips HealthSuite health app on your smartphone so you can monitor changes over time. Compatible with iOS 9.3 or higher/Android 4.4 KitKat or higher. For an overview of (in)compatible devices, please go to [www.philips.com/personalhealthcompatibility](http://www.philips.com/personalhealthcompatibility).

## Supports 2 users



Allows 2 different users to measure their blood pressure levels and heart rate, and monitor changes over time in the Philips health app.

## Rechargeable



No batteries required: rechargeable using USB cable

## Easy-to-read display



Large backlit screen clearly shows your measurement results.

## Easy-to-carry



Easy-to-carry portable wrist unit allows you to monitor your blood pressure and heart rate anywhere

# Specifications

## Design

Colour: Black

Device dimensions: 85 x 67 x 19 mm

Display: LCD display, white backlight

Display dimensions: 47 x 29 mm

Net. Weight: 100 g

## Ease of use

Can be used as stand-alone

Auto sync data

Blood pressure classification

Irregular heartbeat detector

Measure during inflation

LCD display

Supports 2 users

## Accessories

Cuff

USB charger

Wrist blood pressure monitor

## Power

Adapter type: 5 V/1 A USB charger

Battery life: 20 days\*

Battery Type: Rechargeable li-polymer batt.

Charging time: ± 2 hrs

## Service

2 year guarantee

## Syncing

Data transfer technology: Bluetooth®

Phone compatibility: iOS9.3+  
(iPhone5+)/Android4.4+

## Technical

Measurement pressure: 5.3 kPa–30.7 kPa  
(40–230 mmHg)

Heart rate accuracy: ±5% of result on display

Pulse value measurement range: 40–199

heart beats/minute

Rated cuff pressure: 0 kPa–40 kPa (0 mmHg–  
300 mmHg)

Cuff size of the upper arm: 13.5–21.5 cm

Pressure accuracy: 5–40 °C within ±0.4 kPa (3  
mmHg)

Operating environment: Atmospheric  
pressure: 86–106 KPa, Relative humidity: <  
85%RH, Temperature: 5 °C to 40 °C

