

Philips
HipCase

for iPhone

DLA1030H



Carry your iPhone in style

Carry your iPod in style

Stylish and durable, the holster-style HipCase combines sleek protection with convenient one-handed access. Featuring a horizontal belt clip, HipCase provides the utmost in comfort and accessibility. Inner pocket holds an ID card or cash.

Attractive protection

- Modern nylon exterior
- Soft fabric-lined interior

Convenient access

- Access to headphone jack
- Centre notch for one-handed device retrieval
- Sturdy belt clip keeps your device close at hand

Easy organisation

- Inner ID pocket

PHILIPS
sense and simplicity

Specifications

Packaging dimensions

- EAN: 87 12581 43228 7
- Number of products included: 1
- Packaging type: Blister
- Packaging dimensions (W x H x D):
13.97 x 20.32 x 3.81 cm
- Gross weight: 0.118 kg
- Tare weight: 0.071 kg
- Net weight: 0.047 kg

Inner Carton

- GTIN: 1 87 12581 43228 4
- Number of consumer packages: 6

- Inner carton (L x W x H):
24.765 x 20.955 x 13.97 cm
- Gross weight: 0.708 kg
- Tare weight: 0.426 kg
- Net weight: 0.282 kg

Outer Carton

- GTIN: 2 87 12581 43228 1
- Number of consumer packages: 36
- Outer carton (L x W x H): 43.6 x 26.6 x 45.2 cm
- Gross weight: 6.1 kg
- Tare weight: 4.408 kg
- Net weight: 1.692 kg



Issue date 2010-03-11

Version: 1.2.4

12 NC: 8670 000 40889
EAN: 87 12581 43228 7

© 2010 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com