




DLO
HipCase Folio

For iPod
classic



DLA1407D



Carry your iPod in style

in a case with convenient clip

Stylish and durable the folio-style case surrounds your iPod in luxurious leather. Use the sturdy belt clip to keep your device close at hand, and flip open the padded cover to reveal an inner pocket - perfect for an ID or cash.

Attractive protection

- Leather exterior for luxurious look and feel
- Soft fabric-lined interior

Convenient access

- Sturdy belt clip keeps your device close at hand
- Access to all controls and dock connector
- Access to headphone jack while still protected

Discreet Storage

- Small inner pocket perfect for ID or cash

Specifications

iPod compatibility

- Compatible with: iPod classic, iPod classic 120GB, iPod classic 160GB, iPod classic 80GB

Design and finishing

- Color(s): Black
- Materials: Leather

Packaging dimensions

- Packaging type: Box
- Number of products included: 1
- Packaging dimensions (W x H x D): 5.0 x 7.0 x 1.6 inch
- Gross weight: 0.201 lb
- Nett weight: 0.099 lb
- Tare weight: 0.101 lb
- UPC: 6 09585 18753 5

Outer Carton

- Number of consumer packagings: 18
- Outer carton (L x W x H): 16.9 x 11.3 x 7.5 inch
- Gross weight: 3.858 lb
- Nett weight: 1.786 lb
- Tare weight: 2.072 lb
- GTIN: 2 06 09585 18753 9

Inner Carton

- Number of consumer packagings: 3
- Inner carton (L x W x H): 5.4 x 5.4 x 7.0 inch
- Gross weight: 0.639 lb
- Nett weight: 0.298 lb
- Tare weight: 0.342 lb
- GTIN: 1 06 09585 18753 2

Highlights

Leather exterior

This case is made out of attractive and durable leather, so it provides protection for your device -- in style.

Soft interior lining

The soft fabric-lined interior of this case protects the finish of your device while it's cradled inside.

Sturdy belt clip

The sturdy belt clip secures your device to your belt, pocket or bag strap so it's always close at hand.