



Philips  
HipCase

For iPod  
nano G3

DLA28000H



## Carry your iPod in style

in a convenient case with belt clip

Stylish and durable, the folio-style HipCase surrounds your iPod in luxurious leather. Use the sturdy belt clip to keep your iPod close at hand, and flip open the padded cover to reveal an inner pocket -- perfect for an ID card or cash.

### Attractive protection

- Leather exterior for luxurious look and feel
- Soft fabric-lined interior

### Convenient access

- Access to all controls and dock connector
- Access to headphone jack while still protected
- Sturdy belt clip keeps your device close at hand

### Easy organization

- Inner ID pocket

# PHILIPS

# Specifications

## Packaging dimensions

- Packaging dimensions (W x H x D):  
11.43 x 20.32 x 3.81 cm
- Gross weight: 0.071 kg
- Tare weight: 0.04 kg
- Nett weight: 0.031 kg
- EAN: 87 12581 43065 8
- Number of products included: 1
- Packaging type: Blister

## Outer Carton

- Outer carton (L x W x H): 34.5 x 25.6 x 41 cm
- Gross weight: 3.1 kg

- Tare weight: 1.984 kg
- Nett weight: 1.116 kg
- GTIN: 2 87 12581 43065 2
- Number of consumer packagings: 36

## Inner Carton

- Inner carton (L x W x H):  
24.13 x 20.95 x 11.43 cm
- Gross weight: 0.425 kg
- Tare weight: 0.239 kg
- Nett weight: 0.186 kg
- GTIN: 1 87 12581 43065 5
- Number of consumer packagings: 6



Issue date 2014-02-05

Version: 3.0.8

12 NC: 8670 000 40687  
EAN: 87 12581 43065 8

© 2014 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)