



Philips  
HipCase

For iPod  
touch G2

DLA71054



## Carry your iPod in style

in a convenient case with belt clip

Stylish and durable, the folio-style HipCase surrounds your iPod in luxurious leather. Use the sturdy belt clip to keep your iPod close at hand, and flip open the padded cover to reveal an inner pocket -- perfect for an ID card or cash.

### Attractive protection

- Leather exterior for luxurious look and feel
- Soft fabric-lined interior

### Convenient access

- Sturdy belt clip keeps your device close at hand
- Access to all controls and dock connector
- Access to headphone jack while still protected

### Easy organization

- Inner ID pocket

# PHILIPS

# Specifications

## Packaging dimensions

- Packaging dimensions (W x H x D):  
13.97 x 20.32 x 3.18 cm
- Nett weight: 0.053 kg
- Gross weight: 0.101 kg
- Tare weight: 0.048 kg
- EAN: 87 12581 46179 9
- Number of products included: 1
- Packaging type: Blister

## Outer Carton

- Outer carton (L x W x H): 43.8 x 24 x 46.2 cm
- Nett weight: 1.908 kg
- Gross weight: 4.7 kg
- Tare weight: 2.792 kg

- GTIN: 2 87 12581 46179 3
- Number of consumer packagings: 36

## Inner Carton

- Inner carton (L x W x H):  
22.225 x 20.32 x 13.97 cm
- Nett weight: 0.318 kg
- Gross weight: 0.605 kg
- Tare weight: 0.287 kg
- GTIN: 1 87 12581 46179 6
- Number of consumer packagings: 6

## Design

- Color: Black



Issue date 2014-02-05

Version: 1.0.3

12 NC: 8670 000 49177  
EAN: 87 12581 46179 9

© 2014 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)