

Philips
HipCase

for iPod classic

DLA91000H



Carry your iPod in style

in a convenient case with belt clip

Stylish and durable, the folio-style HipCase surrounds your iPod in luxurious leather. Use the sturdy belt clip to keep your iPod close at hand, and flip open the padded cover to reveal an inner pocket -- perfect for an ID card or cash.

Attractive protection

- Leather exterior for luxurious look and feel
- Soft fabric-lined interior

Convenient access

- Access to all controls and dock connector
- Access to headphone jack while still protected
- Sturdy belt clip keeps your device close at hand

Easy organization

- Inner ID pocket

PHILIPS

Specifications

Packaging dimensions

- EAN: 87 12581 43070 2
- Number of products included: 1
- Packaging type: Blister
- Packaging dimensions (W x H x D):
13.97 x 20.32 x 3.18 cm
- Gross weight: 0.101 kg
- Tare weight: 0.048 kg
- Nett weight: 0.053 kg

Inner Carton

- GTIN: 1 87 12581 43070 9
- Number of consumer packagings: 6

- Inner carton (L x W x H):
22.225 x 20.32 x 13.97 cm
- Gross weight: 0.605 kg
- Tare weight: 0.287 kg
- Nett weight: 0.318 kg

Outer Carton

- GTIN: 2 87 12581 43070 6
- Number of consumer packagings: 36
- Outer carton (L x W x H): 43.8 x 24 x 46.2 cm
- Gross weight: 4.7 kg
- Tare weight: 2.792 kg
- Nett weight: 1.908 kg



Issue date 2014-02-05

Version: 1.1.6

12 NC: 8670 000 40693
EAN: 87 12581 43070 2

© 2014 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com