

# PHILIPS

Replacement Kit

Performer

4 x dust bags (s-bag® ULP)

1 x HEPA13 exhaust filter

1 x triple-layer motor filter

FC8060/01



## Replacement kit for Performer Pro/Expert/Ultimate

### Original dust bags and filter replacements

Original starter kit for the Philips vacuum cleaners. The kit contains 4 dust bags, together with the motor and exhaust filters which need to be replaced once a year.

#### Philips Original accessories

- s-bag® Ultra Long Performance original dust bags
- HEPA13 exhaust filter for excellent filtration
- Triple-layer motor filter


# Highlights

## s-bag® Ultra Long Performance



The kit contains 4 x s-bag® Ultra Long Performance original dust bags. This dust bag has an XXL 5 liter capacity and therefore it helps maintain the cleaning performance for longer. Made of premium, multilayer, non-clogging material, the s-bag Ultra Long Performance guarantees excellent suction power and filtration, lasting up to 80% longer than the standard bag.

## HEPA13 exhaust filter



The kit contains 1 x HEPA13 exhaust filter. The HEPA13 filter captures more than 99.95% of fine dust, which other filters simply recirculate back into the air of your home. It traps the fine particles, such as pollens and dust mites, which trigger allergy and asthma symptoms. The filter should be replaced once a year.

## Triple-layer motor filter



The kit contains 1 x triple-layer motor filter. This filter prevents fine dust from reaching the motor and damaging it. The filter should be replaced once a year.

# Specifications

## Contains

s-bag Ultra Long Performance: 4x  
HEPA filter: 1x  
Triple-layer inlet filter: 1x

## Suitable for

Jewel: FC9050 - FC9079  
Performer: FC8680 - FC8682, FC9150 - FC9179  
PerformerPro: FC9180 - FC9199  
Performer Expert: FC8720 - FC8728  
Performer Silent: FC8741, FC8743, FC8745, FC8779 - FC8786  
Performer Ultimate: FC8921 - FC8925, FC8941 - FC8957

