

Active Crystals

PHILIPS

SWAROVSKI

USB MEMORY KEY

GIVE SOUL TO YOUR DATA

for documents, pictures and music

Functionality at its ultimate. Fantasy at its brightest. From now on, a sensational connection between Philips and Swarovski brings technology to new heights of enjoyment. Its name: Active Crystals.

FINEST BEAUTY


- Jet hematite crystals set in Ceralun
- Stylish jewel in beautiful polished stainless steel
- Beautiful pouch to protect your precious USB Memory Key

PERFORMANCE AND CAPACITY

- 1 GB storage capacity for large data files
- Fast data transfer with high-speed USB 2.0
- Share and store up to 250 songs and 1000 photos
- Specially designed and easy to use software included

SAFE AND SECURE

- Protect your valuable personal data with a password



PHILIPS AND SWAROVSKI ACTIVE CRYSTALS

USB MEMORY KEY – FM01SW81/00

SPECIFICATIONS

STORAGE MEDIA

- Built-in Memory Capacity: 1 GB
- Built-in memory type: MLC NAND Flash
- Transfer rate: reads max 15 MB/s; writes max 3 MB/s

CONNECTIVITY

- USB: USB 2.0 high speed

ACCESSORIES

- Storage pouch
- User Manual: English, French, German, Spanish, Portuguese, Italian, Russian, Traditional Chinese, Korean, Japanese
- Warranty booklet

SYSTEM REQUIREMENTS

- PC OS: Windows Vista, XP, 2000, ME, 98SE; Mac OS 9.0 and higher; Linux 2.4.0 and higher
- USB: Free USB port

MISCELLANEOUS

- Guaranteed: 2 year

PACKAGING DATA

- Packaging type: Jewel case in beauty box

OUTER CARTON

- Gross weight: 18.695 lb
- Gross weight: 8.48 kg
- Outer carton (L x W x H): 19.5 x 9.4 x 10.9 inch
- Outer carton (L x W x H): 49.5 x 24 x 27.6 cm
- Net weight: 2.469 lb
- Net weight: 1.12 kg
- Tare weight: 16.226 lb
- Tare weight: 7.36 kg
- EAN: 87 12581 39112 6
- Number of consumer packages: 16

PACKAGING DIMENSIONS

- Packaging dimensions (W x H x D): 3.8 x 3.8 x 3.8 inch
- Packaging dimensions (W x H x D): 9.6 x 9.6 x 9.6 cm
- Gross weight: 0.816 lb
- Gross weight: 0.37 kg

- Net weight: 0.154 lb
- Net weight: 0.07 kg
- Tare weight: 0.661 lb
- Tare weight: 0.3 kg
- EAN: 87 12581 39111 9
- Number of products included: 1
- Packaging type: Carton

INNER CARTON

- Gross weight: 4.674 lb
- Gross weight: 2.12 kg
- Inner carton (L x W x H): 9.4 x 9.0 x 5.0 inch
- Inner carton (L x W x H): 23.9 x 22.8 x 12.6 cm
- Net weight: 0.617 lb
- Net weight: 0.28 kg
- Tare weight: 4.056 lb
- Tare weight: 1.84 kg
- EAN: 87 12581 39113 3
- Number of consumer packages: 4

PRODUCT HIGHLIGHTS

CERALUN JET HEMATITE CRYSTALS

Style and Simplicity gracefully fuse with crystal light, due to the Jet hematite crystals set in Ceralun. When useful is not enough, the new Active Crystals USB Memory Key gives soul to your data sharing.

POLISHED STAINLESS STEEL BODY

Functionality at its ultimate and fantasy at its brightest. The new Active Crystals USB Memory Key comes in a stylish and precious design, made of polished stainless steel.

BEAUTIFUL STORAGE POUCH

Attach your Active Crystals USB Memory Key to your bag and wear it as a stylish charm. To protect and store your precious Active Crystals USB Memory Key, a specially designed pouch is included.

1 GB STORAGE CAPACITY

A useful 1 GB gives you the capacity to store, swap and share large files like documents, music and pictures via the USB port of your PC or laptop.

FAST DATA TRANSFER

High-speed data transfer rates significantly cut annoying waiting time when you are copying large multimedia files to or from your computer's hard disk.

UP TO 250 SONGS AND 1000 PHOTOS


This stylish jewel has been developed to share and store up to 250 songs and 1000 pictures and makes it ideal for personal, portable and powerful use. From a mother carrying around pictures of her children, to business women carrying their files, with extreme ease and convenience.

SOFTWARE INCLUDED

A complete software suite is included, which will run automatically from the Active Crystals USB Memory Key and give you access to useful features, like file synchronisation and automatic file compression.

PASSWORD PROTECTION

Due to password protection including 256-bit AES encryption, your valuable personal data can be protected from unauthorised access by simply using a password.



Issue date 2011-09-23
Version: 4.0.6
12 NC: 8670 000 36526
EAN: 87 12581 39111 9

© 2011 Koninklijke Philips Electronics N.V. and Swarovski AG. All Rights reserved.

Specifications are subject to change without notice.

Trademarks are the property of Koninklijke Philips Electronics N.V., Swarovski AG or their respective owners.

www.active-crystals.com