

CD Mini System

FW775P

Dolby® Pro Logic 3 CD Changer Mini System

- 2 x 240 watts + 80 watts + 2 x 40 watts music power
- 2 x 120 watts + 40 watts + 2 x 20 watts RMS Stereo
- 5 speaker Dolby® Pro Logic Surround Sound System
- Incredible Surround™
- G-Sound Technology
- Dual full logic tape deck
- Auto reverse (R+1)
- Digital sound control (Optimal, Jazz, Rock, Techno, Classic)
- 3-step dynamic bass boost (Beat, Punch, Blast)
- Subwoofer ready
- 40 track CD programmability
- FM/MW/LW stereo digital tuner with 40 presets
- Full RDS suite including "NEWS"
- 4 way bass reflex speaker system with detachable grilles
- Easy set for tuner and clock
- Bedroom clock with wake up alarm
- Wireless prepared

2x240+80+2x40
watts MPO



Let's make things better.



PHILIPS

CD Mini System FW775P

Dolby® Pro Logic 3 CD changer mini system

standard product information

amplifier

- 2 x 240 watts + 80 watts + 2 x 40 watts music power
- 2 x 120 watts + 40 watts + 2 x 20 watts RMS, ±10% THD, 6 ohm, 1 kHz
- 5 Speaker Dolby® Pro Logic Surround Sound System
- Incredible Surround™
- digital sound control (optimal, jazz, rock, techno, classic)
- 3-step dynamic bass boost (beat, punch, blast)
- signal attenuated bass control

compact disc player

- motorized front loading 3 CD changer
- digital servo technology
- bit-check D/A converter
- 40 track CD random program
- play mode (repeat/shuffle/program)
- change 2 discs while play 1

tuner

- RDS - station name, programme type, RDS text, clock set
- RDS - "NEWS" feature
- FM/MW/LW stereo digital tuner with 40 presets
- easy set/auto store for tuner and clock

tape deck

- dual full logic tape deck
- auto reverse (R+1)
- high speed dubbing
- CD synchro recording
- continuous play
- automatic recording level control
- automatic tape type selection
- electronic speed control and full auto stop

optional accessories :



FB201 Active Subwoofer



FB202W Wireless Active Subwoofer

loudspeaker boxes

- 4 way bass reflex
- 6.5" G-cone woofer
- 2.5" gradually cooled ferro fluid tweeter
- 2 x 1" polydome tweeters
- magnetically shielded
- detachable grilles
- hyper energy bass ports for deep powerful bass reproduction

surround speaker

- 1 way speaker (4" full range)
- 6 Ohm

center speaker


- 1 way speaker (4" full range)
- 12 Ohm
- magnetically shielded

connections

- stereo headphones : 3.5mm stereo jack
- video/aux : cinch type sockets
- loudspeaker : click-fit connections
- subwoofer output (see optional accessories)
- line out - wireless prepared (see optional accessories)
- surround : click fit and cinch type

general

- fluorescent tube display
- clock/timer
- wake up with CD, tuner or tape
- 32-key remote control
- PINGS (Philips Intelligent Navigation Guidance System) for ease of operation
- automatic clock set via RDS



PHILIPS

/22 /34

Europe

E Europe

tuner

FM65.81-74MHz

FM87.5-108MHz

MW522-1602kHz

LW153-279kHz

FM75 ohm clickfit

RDS

mains supply

230V/50Hz

general

made in sticker


G-Sound Technology

Philips invented Gyros switching. Controlled music power pumped from the super class G amplifier to speaker enclosures with tuned Bass Ports, gradually cooled ferro-fluid, hybrid polydome tweeters and woofers with G-cone speaker technology.



5 Speaker Surround Sound System

Philips 5 Speaker Dolby® Pro Logic Surround Sound System treats homeviewers to cinema sound quality that has until now, remained in the realm of cinemas. Simply connect to hiFi VCR or stereo TV.



Product Dimensions:

267mmV x 310mmH x 310mmD (electronics)

240mmV x 310mmH x 305mmD (speakers)

200mmV x 98mmH x 98mmD (surround)

260mmV x 140mmH x 240mmD (center)

Package Dimensions:

777mmV x 410mmH x 654mmD

Product weight: 23.5 kg

Carton contains: CD Mini System, Surround Sound Speakers and Center Speaker, Multi-Function Remote Control, Batteries, Operating instructions, Warranty information, FM antenna, AM antenna

Printed in Singapore

Specifications subject to change without notice.

Incredible Surround™ is a Trademark of Philips Consumer Electronics Company

Dolby® Pro Logic is a trademark of Dolby Laboratories

Let's make things better.



PHILIPS