

PHILIPS

NanoProtect filter -
annual filter set

Series 2000

Captures 99.97% of particles
Reduces allergens, gases, odor

FY2423/40



Intelligent purification

Effectively reduces allergens, gases and odors

The Philips NanoProtect annual filter set includes one True HEPA filter and 2 Active Carbon filters. It is compatible with the Philips Air Purifier Series 2000 and 2000i.

Superior purification

- True HEPA filter captures 99.97% particles
- Active carbon filter effectively reduces TVOC and odors

Long-lasting performance

- Replace True HEPA filter approx. every 12 months
- Replace Active Carbon filter approx. every 6 months

Designed for convenience

- 3 filters included - enough for approximately 1 year supply

Highlights

3 Filters included

The Philips NanoProtect annual filter set includes one True HEPA filter and 2 Active Carbon filters. It is enough for approximately 1 year supply.

reduces gases and odors

Philips NanoProtect Active Carbon filter is proven to reduce gases and odors.

Active Carbon lifetime

The Active Carbon filter only needs to be replaced approximately every 6 months.*

True HEPA filter lifetime

The True HEPA filter only needs to be replaced approximately every 12 months.*

99.97% filtration

Philips NanoProtect True HEPA filter captures up to 99.97% of particles that pass through the filter as small as 0.3 microns – the size of most common airborne allergens, particles, bacteria and viruses.

Specifications

Weight and dimensions

Weight of product: 0,879 kg

Weight incl. packaging: 1,134 kg

Dimension of product (LxWxH): 362 x 276 x 70

Dimension of packaging (LxWxH): 79 x 294 x 371 mm

Country of origin

Made in: China

Logistic data

12NC code: 883442340540

