

# PHILIPS

NanoProtect filter -  
annual filter set

Series 2000

Captures 99.97% of particles  
Reduces allergens, gases, odor

FY2423/40



## Intelligent purification

Effectively reduces allergens, gases and odors

The Philips NanoProtect annual filter set includes one True HEPA filter and 2 Active Carbon filters. It is compatible with the Philips Air Purifier Series 2000 and 2000i.

### Superior purification

- True HEPA filter captures 99.97% particles
- Active carbon filter effectively reduces TVOC and odors

### Long-lasting performance

- Replace True HEPA filter approx. every 12 months
- Replace Active Carbon filter approx. every 6 months

### Designed for convenience

- 3 filters included - enough for approximately 1 year supply

# Highlights

## 3 Filters included

The Philips NanoProtect annual filter set includes one True HEPA filter and 2 Active Carbon filters. It is enough for approximately 1 year supply.

## reduces gases and odors

Philips NanoProtect Active Carbon filter is proven to reduce gases and odors.

## Active Carbon lifetime

The Active Carbon filter only needs to be replaced approximately every 6 months.\*

## True HEPA filter lifetime

The True HEPA filter only needs to be replaced approximately every 12 months.\*

## 99.97% filtration

Philips NanoProtect True HEPA filter captures up to 99.97% of particles that pass through the filter as small as 0.3 microns – the size of most common airborne allergens, particles, bacteria and viruses.

# Specifications

## Weight and dimensions

Weight of product: 0,879 kg

Weight incl. packaging: 1,134 kg

Dimension of product (LxWxH): 362 x 276 x 70

Dimension of packaging (LxWxH): 79 x 294 x 371 mm

## Country of origin

Made in: China

## Logistic data

12NC code: 883442340540

