

# PHILIPS

## Anti-calc cartridges

GC018



## Anti-calc cartridges

### Extend the lifetime of your Wardrobe Care

The Wardrobe Care is equipped with a unique DualProtect Anti-calc system, which keeps the appliance 99% calc-free. The special anti-calc cartridges demineralise the water. You can enjoy ironing for longer, without having to worry about calc.

#### **Enhance the lifetime of your Wardrobe Care**

- Anti-calc cartridges
- Keeps your ironing system 99% calc-free
- Suitable for WardrobeCare with Optimal TEMP (GC9940, GC9950)

#### **Easy to use**

- Easy to replace
- Replace cartridge once the granules have changed colour
- Easily store spare cartridges in the special storage compartment
- Two anti-calc cartridges in one pack

# Highlights


## Easy to replace

The anti-calc cartridge is easy to replace. Simply take out the anti-calc cartridge from the water tank and put in a new one. The old anti-calc cartridge can be disposed off with household rubbish.

## Anti-calc cartridges

Anti-calc cartridges demineralise the water and prevent calc from entering the appliance.

## Easy to store cartridges



The Wardrobe Care integrated ironing solution has a special storage compartment where spare anti-calc cartridges can easily be stored.

## Keep your system 99% calc-free



The anti-calc cartridges keep your Wardrobe Care free from calc.

## Replace cartridges

The anti-calc cartridges will need to be replaced every 3 to 7 months. Once the granules in the cartridge have completely changed colour, you will need to replace it.

## Suitable for GC9940, GC9950

These anti-calc cartridges are suitable for Wardrobe Care integrated ironing solutions with type numbers GC9940 and GC9950.

## 2 cartridges per pack

One package contains two anti-calc cartridges.

# Specifications

## Calc management

Suitable for tap water

## Accessories

Suitable for Wardrobe Care: GC9940, GC9950

